

FIG. 1

102




FIG. 2

Unit Operations List

Microbial Fermentation Process

UOP Seq. No.	Code	Unit Operation Type	Cycles per		Batch				Process				Recovery			
			UnOp Offset (Hrs)		UnOp Start	UnOp End	Offset (Hrs)		UnOp Start	UnOp End	Offset (Hrs)		SWR	Product	OAR	Total Protein
1	1	Inoculum Prep	1	3	1	6			1							
2	2	Flask Growth	1	3	1	6			1							
3	53	Seed Fermentation	1	3	1	6			1							
4	3	Production Fermentation	1	3	1	6			1							
5	51	Heat Exchange	1	3	1	6			1							
6	28	Cont. Centrifugation/Whole Cell Harvest	1	3	1	6			1							
7	48	Resuspend Cell Paste	1	1	1	6			1							
8	51	Heat Exchange	1	3	8	10			1							
9	31	Cell Disruption/ High Pressure	1	3	8	10			1							
10	51	Heat Exchange	1	3	8	10			1							
11	48	Resuspension/ Surfactant	1	2	11	12			1							
12	29	Cont. Centrifugation/ Precipitate Harvest	1	2	11	12			1							
13	48	Resuspension/ Buffer	1	1	1	12			1							
14	29	Ultrafiltration/ Concentration/ Dilution	1	1	1				1							
15	48	Microfiltration/ Tangential Flow	1	1	1				1							
16	36	Product Adsorption MPLC	1	1	1				1							
17	34	Product Adsorption MPLC	1	1	1				1							
18	39	Ultrafiltration/ Flow Dialysis	1	1	1				1							
19	39	Product Adsorption MPLC	1	1	1				1							
20	37	Ultrafiltration/ Flow Dialysis	1	1	1				1							
21	39	Product Adsorption MPLC	1	1	1				1							
22	37	Microfiltration/ Dead End	1	1	1				1							
23	99	End	1	1	1				1							
			306		312		314		316		318		320		322	
			308		310		312		314		316		318		320	
			304		306		308		310		312		314		316	
			302		304		306		308		310		312		314	
			328		330		332		334		336		338		340	

FIG. 3

1. The first part of the document is a list of references. The references are listed in a standard format, with the author's name, the title of the work, and the publisher. The references are as follows:

Mammalian Cell Culture Process

UOP Seq. No.	Code	Unit Operation Type	Cycles per			Batch			Process			Recovery		
			UnOp		Offset (Hrs)	UnOp		UnOp End	UnOp End	UnOp End	Product		Total Protein	
			Offset (Hrs)	UnOp Start		UnOp End	SWR				OAR	SWR	OAR	
1	4	Initial Seeding	1			1								
2	5	Culture Vessel Split	1			1								
3	5	Culture Vessel Split	1			1								
4	5	Culture Vessel Split	1			1								
5	6	Spinner Flask Split	1			1								
6	54	Spinner Flask Split	1			1								
7	13	Stirred Tank Reactor	1			1								
8	61	Harvest/Feed	7	24		1			8	18				168
9	62	Harvest Pool	1			1			8	18				168
10	34	MF/Tangential Flow	1			1			8	18				168
11	36	UF/Concentration	1			1			8	18				168
12	39	PAC/MPLC	1			1			8	18				168
13	39	PAC/MPLC	1			1			8	18				168
14	36	UF/Concentration	1			1			8	18				168
16	39	PAC/MPLC	1			1			8	18				168
16	37	UF/Flow Dialysis	1			1			8	18				168
17	39	PAC/MPLC	1			1			8	18				168
18	35	MF/Dead End	1			1			8	18				168
19	99	End	1			1			8	18				168

406 408 410 412 414 416 418 420 422 424

402 404

FIG. 4




FIG. 5

unit operation id code	Unit operation type	Parameters	solution type	tasks	task duration
1	inoculum prep	# of flasks, volume of flasks, temperature, agitation, duration, final OD	S-101	setup, preincubation, incubation, clean up	3, 3, 23, .3 Hrs
2	flask growth	scale up ratio, media volume, temperature, agitation, duration, final OD	S-101	setup, preincubation, incubation, clean up	1, 1, 23, .3 Hrs
3	fermentation seed	scale up ratio, fermentor working volume, antifoam, base acid, grow temperature, agitation, sparge rate, back pressure, total duration	S-101, 102, 103, 104, 105	setup, preincubation, fermentation, harvest, CIP, SIP, clean up	1, 1, 21, .5, 1, 1, 3 Hrs
4	fermentation production	scale up ratio, fermentor working volume, antifoam A, antifoam B, base, acid, grow temperature, agitation, sparge rate, back pressure, total duration, final OK, dry cell mass, product concentration, CIP, SIP	S-101, 102, 103, 104, 105	setup preincubation, fermentation, CIP, SIP, cleanup	.
5	heat exchange	process initial & final temp; utility initial & final temp; process specific heat; design type, step recovery of product, step recovery of T.P., temperature regulation, CIP, SIP		setup, transfer, CIP, SIP, cleanup	.
6	batch centrifugation	system void volume, RCF, time, volume reduction, wash volume, clean, rinse	S-105	setup, centrifugation, wash, CIP, SIP, cleanup	.
7	resolubilization resuspension	reagent/product ratio, titration solution, resolubilization, agitation, solution name, step recovery of the product, step recovery of T.P., temperature regulation, CIP, SIP	S-107	setup, dilution, agitate, CIP, SIP, clean up	.
8	Cell Disruption High Press. Homogenization	product temperature, utility temperature, void volume, number of passes, pressure, flow rate, temperature increase, wash, rinse, step recovery of product, step recovery of T.P., temperature regulation, CIP	S-107	setup, lysis, CIP, SIP, clean up	.
9	Dilute with Surfactant	reagent product ratio, titration solution, dilution time, agitation, solution name, step recovery of product, step recovery of T.P., temperature regulation, CIP, SIP	S-108	setup, dilution, agitate, CIP, SIP, clean up	.
10	batch centrifugation precipitate harvest	system void volume, RCF, time, volume reduction, wash volume, clean, rinse, step recovery of product, step recovery of T.P., temperature regulation, CIP, SIP	S-108	setup, centrifugation, wash, CIP, SIP, clean up	.
11	resuspend with chaotrope	reagent/product ratio, titration solution, resolubilization, agitation, solution name, step recovery of product, step recovery to TP, temperature regulation, CIP, SIP	S-109	setup, flush, prime, concentration, dilution, wash, flush, store, CIP, SIP, cleanup	.
.
.
.
.

Fig. 6

108
=

2025-03-08 10:00:00




FIG. 7




FIG. 9

Sample Application of Process Design Cycles In Process Scheduling

Microbial Fermentation Process (see unit operation list)

Duration	First Process Cycle		Second Process Cycle	
	Week	Day	Week	Day

Note: None of the unit operations in this process have more than 1 cycle per unit operation
(see unit operation 8 in the mammalian cell culture process for an example of multiple cycles per unit operation)

Unit Operations 1-6 undergo three repetitive cycles per batch as a set before continuing with unit op 7
This translates to three runs on a fermentor with each harvest (unit op 5 & 6) being stored for pooling at unit op 7
Associated with each fermentor run (unit op 4) are the previous steps for inoculation prep (unit ops 1-3)

1/3 fermentation cycles per batch

1	Inoculum Prep	24 hrs	1	Fri - Sat	2	Fri - Sat
2	Flask Growth	24 hrs	2	Sat - Sun	3	Sat - Sun
3	Seed Fermentation	24 hrs	2	Sun - Mon	3	Sun - Mon
4	Production Fermentation	24 hrs	2	Mon - Tue	3	Mon - Tue
5	Heat Exchange	1 hr	2	Tue	3	Tue
6	Centrifugation	1hr	2	Tue	3	Tue

2/3 fermentation cycles per batch

1	Inoculum Prep	24 hrs	2	Sun - Mon	3	Sun - Mon
2	Flask Growth	24 hrs	2	Mon - Tue	3	Mon - Tue
3	Seed Fermentation	24 hrs	2	Tue - Wed	3	Tue - Wed
4	Production Fermentation	24 hrs	2	Wed - Thu	3	Wed - Thu
5	Heat Exchange	1 hr	2	Thu	3	Thu
6	Centrifugation	1hr	2	Thu	3	Thu

3/3 fermentation cycles per batch

1	Inoculum Prep	24 hrs	2	Tue - Wed	3	Tue - Wed
2	Flask Growth	24 hrs	2	Wed - Thu	3	Wed - Thu
3	Seed Fermentation	24 hrs	2	Thu - Fri	3	Thu - Fri
4	Production Fermentation	24 hrs	2	Fri - Sat	3	Fri - Sat
5	Heat Exchange	1 hr	2	Sat	3	Sat
6	Centrifugation	1hr	2	Sat	3	Sat

Unit Operation 7 pools the harvests from the three fermentation cycles above

7	Pool Harvests	3 hr	3	Mon	4	Mon
---	---------------	------	---	-----	---	-----

Unit Operations 8-9 undergo three repetitive cycles per batch as set before continuing with unit operation 11
This translates to three consecutive passes through cell disruptor (unit op 9) with its associated heat exchangers
(unit op 8 & 10) at the inlet and the outlet of the cell disruptor

1/3 disruption cycles per batch

8	Heat Exchange					
9	Cell Disruption					
10	Heat Exchange	0.5 hr	3	Mon	4	Mon

2/3 disruption cycles per batch

8	Heat Exchange					
9	Cell Disruption					
10	Heat Exchange	0.5 hr	3	Mon	4	Mon

3/3 disruption cycles per batch

8	Heat Exchange					
9	Cell Disruption					
10	Heat Exchange	0.5 hr	3	Mon	4	Mon

Sample Application of Process Design Cycles in Process Scheduling

Microbial Fermentation Process (see unit operation list)

			First Process Cycle		Second Process Cycle	
		Duration	Week	Day	Week	Day
Unit ops 11-12 undergo two repetitive cycles per batch as a set before continuing with unit op 13 This translates to two cycles of resuspending the cell lysate from the cell disruptor in a mild surfactant and reconcentrating the insoluble product to a paste by centrifugation						
1/2 product washing cycles per batch						
11	Resuspension	0.5 hr	3	Mon	4	Mon
12	Centrifugation	1 hr	3	Mon	4	Mon
2/3 product washing cycles per batch						
11	Resuspension	0.5 hr	3	Mon	4	Mon
12	Centrifugation	1 hr	3	Mon	4	Mon
Unit ops 13-22 undergo only one cycle per unit operation each to the end of the process						
13	Resuspension	0.5 hr	3	Mon	4	Mon
14	Buffer Exchange	2 hr	3	Mon	4	Mon
15	Filtration	2 hr	3	Mon	4	Mon
16	Liquid Chromatography	16 hrs	3	Mon - Tue	4	Mon - Tue
17	Liquid Chromatography	4 hrs	3	Tue	4	Tue
18	Buffer Exchange	2 hrs	3	Tue	4	Tue
19	Liquid Chromatography	2 hrs	3	Wed	4	Wed
20	Buffer Exchange	2 hrs	3	Wed	4	Wed
21	Liquid Chromatography	2 hrs	3	Wed	4	Wed
22	Filtration	2 hrs	3	Wed	4	Wed

Operation		Duration (Hrs)		Rel. Time Scale (Hrs)		Abs. Days		Start		Finish		Calculations
		Calc.	Adj.	Prep	Exec.	Compil.	Start	End	Date	Date	Time	
1 A Inoculum Prep					15.5		06/03/96	08:00 AM				
1	Set Up	3.0	0.0	3.0 Hrs								
2	Preincubation	3.0	0.0	3.0 Hrs	12.5		06/03/96	09:30 AM	06/03/96	12:30 PM		
3	Incubation	23.0	0.0	23.0 Hrs	15.5		06/03/96	12:30 PM	06/03/96	02:30 PM		
4	Clean Up	0.3	0.0	0.3 Hrs	38.5		06/03/96	02:30 PM	06/03/96	02:30 PM		
5	Subtotal	29.0		29.0 Hrs	38.5		06/03/96	02:30 PM	06/03/96	02:30 PM		
2 A Flask Growth												
6	Set Up	1.0	0.0	1.0 Hrs	37.5		06/04/96	12:30 PM	06/04/96	01:30 PM		
7	Preincubation	1.0	0.0	1.0 Hrs	38.5		06/04/96	01:30 PM	06/04/96	02:30 PM		
8	Incubation	23.0	0.0	23.0 Hrs	61.5		06/04/96	02:30 PM	06/04/96	01:30 PM		
9	Clean Up	0.3	0.0	0.3 Hrs	61.8		06/04/96	01:30 PM	06/04/96	01:45 PM		
10	Subtotal	25.0		25.0 Hrs	61.5		06/04/96	01:30 PM	06/04/96	01:45 PM		
3 A Seed Fermentation												
11	Set Up	1.0	0.0	1.0 Hrs	60.5		06/05/96	11:30 AM	06/05/96	12:30 PM		
12	Preincubation	1.0	0.0	1.0 Hrs	61.5		06/05/96	12:30 PM	06/05/96	01:30 PM		
13	Fermentation	21.0	0.0	21.0 Hrs	82.5		06/05/96	01:30 PM	06/05/96	10:30 AM		
14	Harvest	0.5	0.0	0.5 Hrs	83.0		06/05/96	10:30 AM	06/05/96	11:30 AM		
15	CIP	1.0	0.0	1.0 Hrs	83.5		06/05/96	11:30 AM	06/05/96	12:30 PM		
16	SIP	1.0	0.0	1.0 Hrs	84.5		06/05/96	12:30 PM	06/05/96	01:30 PM		
17	Clean Up	3.0	0.0	3.0 Hrs	87.5		06/05/96	12:30 PM	06/05/96	03:30 PM		
18	Subtotal	28.5		28.5 Hrs	83.0		06/05/96	12:30 PM	06/05/96	03:30 PM		
4 A Production Fermentation												
19	Set Up	1.0	0.0	1.0 Hrs	82.0		06/06/96	09:00 AM	06/06/96	10:00 AM		
20	Preincubation	1.0	0.0	1.0 Hrs	83.0		06/06/96	10:00 AM	06/06/96	11:00 AM		
21	Fermentation	21.0	0.0	21.0 Hrs	104.0		06/06/96	11:00 AM	06/06/96	08:00 AM		
22	CIP	1.0	0.0	1.0 Hrs	105.0		06/06/96	08:00 AM	06/06/96	09:00 AM		
23	SIP	1.0	0.0	1.0 Hrs	106.0		06/06/96	09:00 AM	06/06/96	10:00 AM		
24	Clean Up	2.0	0.0	2.0 Hrs	108.0		06/06/96	10:00 AM	06/06/96	12:00 PM		
25	Subtotal	27.0		27.0 Hrs	104.0		06/06/96	10:00 AM	06/06/96	12:00 PM		
5 A Heat Exchange												
26	Set Up	0.50	0.0	0.5 Hrs	104.5		06/07/96	08:00 AM	06/07/96	08:30 AM		
27	Transfer	1.00	0.0	1.0 Hrs	105.0		06/07/96	09:00 AM	06/07/96	10:00 AM		
28	CIP	1.0	0.0	1.0 Hrs	106.0		06/07/96	10:00 AM	06/07/96	11:00 AM		
29	SIP	1.0	0.0	1.0 Hrs	107.0		06/07/96	11:00 AM	06/07/96	12:00 PM		
30	Clean Up	2.0	0.0	2.0 Hrs	109.0		06/07/96	12:00 PM	06/07/96	01:00 PM		
31	Subtotal	5.0		5.0 Hrs	105.0		06/07/96	11:00 AM	06/07/96	01:00 PM		
6 A Cont. Cant./Solids												
32	Set Up	1.00	0.0	1.0 Hrs	105.0		06/07/96	08:00 AM	06/07/96	09:00 AM		
33	Centrifugation	1.00	0.0	1.0 Hrs	106.0		06/07/96	09:00 AM	06/07/96	10:00 AM		
34	Wash	0.10	0.0	0.1 Hrs	106.1		06/07/96	10:00 AM	06/07/96	10:30 AM		
35	CIP	0.25	0.0	0.3 Hrs	108.4		06/07/96	10:30 AM	06/07/96	11:00 AM		
36	SIP	1.00	0.0	1.0 Hrs	107.4		06/07/96	11:00 AM	06/07/96	12:00 PM		
37	Clean Up	0.50	0.0	0.5 Hrs	107.9		06/07/96	12:00 PM	06/07/96	12:30 PM		
38	Subtotal	3.85		3.85 Hrs	108.1		06/07/96	11:21 AM	06/07/96	11:51 AM		
1 B Inoculum Prep												
39	Set Up	1.0	0.0	1.0 Hrs	14.5		06/03/96	01:30 PM	06/03/96	02:30 PM		
40	Preincubation	1.0	0.0	1.0 Hrs	15.5		06/03/96	02:30 PM	06/03/96	03:30 PM		

15

14

13

12

11

10

9

8

7

6

5

4

3

2

1

FIG-12A

Process Time Line										Rel. Time Scale (Hrs)				Abs. Days				Start				Finish				Calculations			
Duration (Hrs.)										Calc.		AD		Adj.		Prep		Exec.		Compl.		Start		Date		Time		Time	
Operation																													
59	Incubation									23.0	0.0	0.0	23.0	Hrs				15.5				0.65	1.60	06/03/98	03:30 AM	06/04/98	02:30 PM		
60	Clean Up									0.3	0.0	0.0	0.3	Hrs				38.5				38.8	1.61	06/04/98	02:30 PM	06/04/98	02:45 PM		
61	Subtotal									25.0			25.0	Hrs				38.5											
2 B Flask Growth																													
62	Set Up									1.0	0.0	0.0	1.0	Hrs				37.5				1.52	1.58	06/04/98	12:30 PM	06/04/98	01:30 PM		
63	Preincubation									1.0	0.0	0.0	1.0	Hrs				38.5				1.56	1.60	06/04/98	01:30 PM	06/04/98	02:30 PM		
64	Incubation									23.0	0.0	0.0	23.0	Hrs				61.5				1.60	2.68	06/04/98	02:30 PM	06/05/98	01:30 PM		
65	Clean Up									0.3	0.0	0.0	0.3	Hrs							61.8	2.57	06/05/98	01:30 PM	06/05/98	01:45 PM			
66	Subtotal									25.0			25.0	Hrs				61.5											
3 B Seed Fermentation																													
67	Set Up									1.0	0.0	0.0	1.0	Hrs				60.5				2.48	2.52	06/05/98	11:30 AM	06/05/98	12:30 PM		
68	Preincubation									1.0	0.0	0.0	1.0	Hrs				61.5				2.52	2.56	06/05/98	12:30 PM	06/05/98	01:30 PM		
69	Fermentation									21.0	0.0	0.0	21.0	Hrs				82.5				2.56	3.44	06/05/98	01:30 PM	06/06/98	10:30 AM		
70	Harvest									0.5	0.0	0.0	0.5	Hrs				83.0				3.44	3.48	06/06/98	10:30 AM	06/06/98	11:00 AM		
71	CIP									1.0	0.0	0.0	1.0	Hrs							83.5	3.44	06/06/98	10:30 AM	06/06/98	11:30 AM			
72	SIP									1.0	0.0	0.0	1.0	Hrs							84.5	3.48	06/06/98	11:30 AM	06/06/98	12:30 PM			
73	Clean Up									3.0	0.0	0.0	3.0	Hrs				83.0				87.5	3.52	06/06/98	12:30 PM	06/06/98	03:30 PM		
74	Subtotal									28.5			28.5	Hrs				83.0											
4 B Production Fermentation																													
75	Set Up									1.0	0.0	0.0	1.0	Hrs				82.0				3.38	3.42	06/06/98	09:00 AM	06/06/98	10:00 AM		
76	Preincubation									1.0	0.0	0.0	1.0	Hrs				83.0				3.42	3.46	06/06/98	10:00 AM	06/06/98	11:00 AM		
77	Fermentation									21.0	0.0	0.0	21.0	Hrs				104.0				3.46	4.33	06/06/98	11:00 AM	06/07/98	08:00 AM		
78	CIP									1.0	0.0	0.0	1.0	Hrs							105.0	4.33	06/07/98	08:00 AM	06/07/98	09:00 AM			
79	SIP									1.0	0.0	0.0	1.0	Hrs							106.0	4.38	06/07/98	09:00 AM	06/07/98	10:00 AM			
80	Clean Up									2.0	0.0	0.0	2.0	Hrs				104.0				106.0	4.42	06/07/98	10:00 AM	06/07/98	11:00 AM		
81	Subtotal									27.0			27.0	Hrs				104.0											
5 B Heat Exchange																													
82	Set Up									0.50	0.0	0.0	0.5	Hrs				104.5				4.33	4.35	06/07/98	08:30 AM	06/07/98	09:30 AM		
83	Transfer									1.00	0.0	0.0	1.0	Hrs				105.0				4.33	4.38	06/07/98	09:00 AM	06/07/98	09:30 AM		
84	CIP									1.0	0.0	0.0	1.0	Hrs							106.0	4.38	06/07/98	09:30 AM	06/07/98	10:00 AM			
85	SIP									1.0	0.0	0.0	1.0	Hrs							107.0	4.42	06/07/98	10:00 AM	06/07/98	11:00 AM			
86	Clean Up									2.0	0.0	0.0	2.0	Hrs				105.0				107.0	4.46	06/07/98	10:30 AM	06/07/98	11:00 PM		
87	Subtotal									5.0			5.0	Hrs				105.0											
6 B Cont. Cent/Solids																													
88	Set Up									1.00	0.0	0.0	1.0	Hrs				105.0				4.33	4.38	06/07/98	08:30 AM	06/07/98	09:30 AM		
89	Centrifugation									1.00	0.0	0.0	1.0	Hrs				106.0				4.38	4.42	06/07/98	09:00 AM	06/07/98	10:00 AM		
90	Wash									0.10	0.0	0.0	0.1	Hrs				106.1				4.42	4.42	06/07/98	09:30 AM	06/07/98	10:06 AM		
91	CIP									0.25	0.0	0.0	0.3	Hrs							106.4	4.42	06/07/98	10:06 AM	06/07/98	10:21 AM			
92	SIP									1.00	0.0	0.0	1.0	Hrs							107.4	4.43	06/07/98	10:21 AM	06/07/98	11:21 AM			
93	Clean Up									0.50	0.0	0.0	0.5	Hrs				106.1				107.9	4.47	06/07/98	11:21 AM	06/07/98	11:51 AM		
94	Sub Total									3.85			3.85	Hrs				106.1											
1 C Inoculum Prep																													
95	Set Up									1.0	0.0	0.0	1.0	Hrs				14.5				0.66	0.60	06/03/98	01:30 PM	06/03/98	02:30 PM		
96	Preincubation									1.0	0.0	0.0	1.0	Hrs				15.5				0.60	0.65	06/03/98	02:30 PM	06/03/98	03:30 PM		
97	Incubation									23.0	0.0	0.0	23.0	Hrs				38.5				0.65	1.60	06/03/98	03:30 PM	06/04/98	02:30 PM		
98	Clean Up									0.3	0.0	0.0	0.3	Hrs				38.5				38.8	1.61	06/04/98	02:30 PM	06/04/98	02:45 PM		
99	Subtotal									25.0			25.0	Hrs				38.5											

Fig. 12B

Operation	Duration (Hrs)			Rel. Time Scale (Hrs)			Abs. Days		Start		Finish		Calculations
	Calc.	AD	Adj.	Prep	Exec.	Comp.	Start	End	Date	Time	Date	Time	
2 C Flask Growth					15.5				08/03/98	08:00 AM			
116 Set Up	1.0	0.0	1.0 Hrs	37.5			1.52	1.58	08/04/98	12:30 PM	08/04/98	01:30 PM	
117 Preincubation	1.0	0.0	1.0 Hrs	38.5			1.58	1.60	08/04/98	01:30 PM	08/04/98	02:30 PM	
118 Incubation	23.0	0.0	23.0 Hrs		61.5		1.60	2.58	08/04/98	02:30 PM	08/05/98	01:30 PM	
119 Clean Up	0.3	0.0	0.3 Hrs			61.5	2.58	2.57	08/05/98	01:30 PM	08/05/98	01:45 PM	
120 Subtotal	25.0		25.0 Hrs										
3 C Seed Fermentation													
121 Set Up	1.0	0.0	1.0 Hrs	60.5			2.48	2.52	08/05/98	11:30 AM	08/05/98	12:30 PM	
122 Preincubation	1.0	0.0	1.0 Hrs	61.5			2.52	2.58	08/05/98	12:30 PM	08/05/98	01:30 PM	
123 Fermentation	21.0	0.0	21.0 Hrs		82.5		2.58	3.44	08/05/98	01:30 PM	08/06/98	10:30 AM	
124 Harvest	0.5	0.0	0.5 Hrs		83.0		3.44	3.48	08/06/98	10:30 AM	08/06/98	11:00 AM	
125 CIP	1.0	0.0	1.0 Hrs			83.5	3.44	3.48	08/06/98	10:30 AM	08/06/98	11:30 AM	
126 SIP	1.0	0.0	1.0 Hrs			84.5	3.48	3.52	08/06/98	11:30 AM	08/06/98	12:30 PM	
127 Clean Up	3.0	0.0	3.0 Hrs			87.5	3.52	3.65	08/06/98	12:30 PM	08/06/98	03:30 PM	
128 Subtotal	28.5		28.5 Hrs		83.0								60.0 L 1.7 LPM = 0.50 Hrs
4 C Production Fermentation													
129 Set Up	1.0	0.0	1.0 Hrs	82.0			3.38	3.42	08/06/98	09:00 AM	08/06/98	10:00 AM	
130 Preincubation	1.0	0.0	1.0 Hrs	83.0			3.42	3.46	08/06/98	10:00 AM	08/06/98	11:00 AM	
131 Fermentation	21.0	0.0	21.0 Hrs		104.0		3.46	4.33	08/06/98	11:00 AM	08/07/98	08:00 AM	
132 CIP	1.0	0.0	1.0 Hrs			105.0	4.33	4.38	08/07/98	08:00 AM	08/07/98	09:00 AM	
133 SIP	1.0	0.0	1.0 Hrs			106.0	4.38	4.42	08/07/98	09:00 AM	08/07/98	10:00 AM	
134 Clean Up	2.0	0.0	2.0 Hrs			108.0	4.42	4.50	08/07/98	10:00 AM	08/07/98	11:00 AM	
135 Subtotal	27.0		27.0 Hrs		104.0		4.42	4.50	08/07/98	10:00 AM	08/07/98	12:00 PM	
5 C Heat Exchange													
136 Set Up	0.50	0.0	0.5 Hrs	104.5			4.33	4.35	08/07/98	08:00 AM	08/07/98	08:30 AM	
137 Transfer	1.00	0.0	1.0 Hrs				4.33	4.38	08/07/98	08:30 AM	08/07/98	09:00 AM	
138 CIP	1.0	0.0	1.0 Hrs			106.0	4.38	4.42	08/07/98	09:00 AM	08/07/98	10:00 AM	
139 SIP	1.0	0.0	1.0 Hrs			107.0	4.42	4.46	08/07/98	10:00 AM	08/07/98	11:00 AM	
140 Clean Up	2.0	0.0	2.0 Hrs			109.0	4.46	4.54	08/07/98	11:00 AM	08/07/98	01:00 PM	
141 Subtotal	5.0		5.0 Hrs		105.0								582.1 L @ 9.4 LPM = 1.00 Hrs
6 C Cont. Cent./Solids													
142 Set Up	1.00	0.0	1.0 Hrs	105.0			4.33	4.38	08/07/98	08:00 AM	08/07/98	09:00 AM	
143 Centrifugation	1.00	0.0	1.0 Hrs				4.38	4.42	08/07/98	09:00 AM	08/07/98	10:00 AM	
144 Wash	0.10	0.0	0.1 Hrs			106.1	4.42	4.43	08/07/98	10:00 AM	08/07/98	10:06 AM	
145 CIP	0.25	0.0	0.3 Hrs			106.4	4.42	4.43	08/07/98	10:06 AM	08/07/98	10:21 AM	
146 SIP	1.00	0.0	1.0 Hrs			107.4	4.43	4.47	08/07/98	10:21 AM	08/07/98	11:21 AM	
147 Clean Up	0.50	0.0	0.5 Hrs			107.8	4.47	4.49	08/07/98	11:21 AM	08/07/98	11:31 AM	
148 Sub Total	3.85		3.85 Hrs		106.1								582.1 L @ 9.4 LPM = 1.00 Hrs 1.0 L @ 0.2 LPM = 0.10 Hrs 20.0 L @ 1.3 LPM = 0.25 Hrs
7 A Recirculation													
149 Set Up	1.00	0.0	1.0 Hrs	106.1			4.38	4.42	08/07/98	09:00 AM	08/07/98	10:00 AM	
150 Dilution	0.50	0.0	0.5 Hrs				4.42	4.44	08/07/98	10:00 AM	08/07/98	10:38 AM	
151 Agle	1.00	0.0	1.0 Hrs			107.8	4.44	4.48	08/07/98	10:38 AM	08/07/98	11:38 AM	
152 CIP	1.00	0.0	1.0 Hrs			108.8	4.48	4.53	08/07/98	11:38 AM	08/07/98	12:38 PM	
153 SIP	1.00	0.0	1.0 Hrs			109.8	4.53	4.57	08/07/98	12:38 PM	08/07/98	01:38 PM	
154 Clean Up	1.00	0.0	1.0 Hrs			110.8	4.57	4.61	08/07/98	01:38 PM	08/07/98	02:38 PM	
155 Subtotal	5.50		5.50 Hrs		107.8								48.7 L @ 1.8 LPM = 0.50 Hrs 0.50 Hrs
8 A Heat Exchange													
156 Set Up	1.00	0.0	1.0 Hrs										
157 Transfer	1.00	0.0	1.0 Hrs										
158 CIP	1.00	0.0	1.0 Hrs										
159 SIP	1.00	0.0	1.0 Hrs										
160 Clean Up	1.00	0.0	1.0 Hrs										
161 Subtotal	5.50		5.50 Hrs										
162													
163													
164													
165													
166													
167													
168													
169													
170													
171													
172													
173													
174													

FIG. 12C

Process Time Line																			
Operation		Duration (Hrs.)		Rel. Time Scale (Hrs)		Abs. Day		Start		Finish		Calculations							
		Calc.	Adj.	Prep	Exec.	Comp.	End	Start	End	Start	End			Time	Time	Time	Time	Time	Time
176	Set Up	0.50	0.0	0.5 Hrs	107.8	15.5	4.46	0603:06	08:00 AM	0607:06	08:00 AM	66.5 L @	3.7 LPM	11:38 AM	11:38 AM				
177	Transfer	0.30	0.0	0.3 Hrs	107.9		4.48	0607:06	11:38 AM	0607:06	11:38 AM			11:54 AM	11:54 AM				
178	CIP	0.0	0.0	0.0 Hrs			4.50	0607:06	11:54 AM	0607:06	11:54 AM			11:54 AM	11:54 AM				
179	SIP	0.0	0.0	0.0 Hrs			4.50	0607:06	11:54 AM	0607:06	11:54 AM			11:54 AM	11:54 AM				
180	Clean Up	0.0	0.0	0.0 Hrs			4.50	0607:06	11:54 AM	0607:06	11:54 AM			11:54 AM	11:54 AM				
181	Subtotal	0.8			107.9														
182	9 A Homogenization																		
183	Set Up	0.25	0.0	0.3 Hrs	107.9		4.49	0607:06	11:39 AM	0607:06	11:39 AM	66.5 L @	1.6 LPM	11:54 AM	11:54 AM			0.68 Hrs	
184	Lysis	0.68	0.0	0.7 Hrs			4.50	0607:06	11:54 AM	0607:06	11:54 AM			12:34 PM	12:34 PM				
185	CIP	0.0	0.0	0.0 Hrs			4.52	0607:06	12:34 PM	0607:06	12:34 PM			12:34 PM	12:34 PM				
186	SIP	0.0	0.0	0.0 Hrs			4.52	0607:06	12:34 PM	0607:06	12:34 PM			12:34 PM	12:34 PM				
187	Clean Up	0.0	0.0	0.0 Hrs			4.52	0607:06	12:34 PM	0607:06	12:34 PM			12:34 PM	12:34 PM				
188	Subtotal	0.9			108.8		4.52	0607:06	12:34 PM	0607:06	12:34 PM								
189	10 A Heat Exchange																		
190	Set Up	0.50	0.0	0.5 Hrs	108.8		4.50	0607:06	12:04 PM	0607:06	12:04 PM	69.0 L @	3.8 LPM	12:34 PM	12:34 PM			0.30 Hrs	
191	Transfer	0.30	0.0	0.3 Hrs			4.52	0607:06	12:34 PM	0607:06	12:34 PM			12:52 PM	12:52 PM				
192	CIP	0.0	0.0	0.0 Hrs			4.54	0607:06	12:52 PM	0607:06	12:52 PM			12:52 PM	12:52 PM				
193	SIP	0.0	0.0	0.0 Hrs			4.54	0607:06	12:52 PM	0607:06	12:52 PM			12:52 PM	12:52 PM				
194	Clean Up	0.0	0.0	0.0 Hrs			4.54	0607:06	12:52 PM	0607:06	12:52 PM			12:52 PM	12:52 PM				
195	Subtotal	0.8			108.9		4.54	0607:06	12:52 PM	0607:06	12:52 PM								
196	8 B Heat Exchange																		
197	Set Up	0.00	0.0	0.0 Hrs	108.9		4.54	0607:06	12:52 PM	0607:06	12:52 PM	66.5 L @	3.7 LPM	01:10 PM	01:10 PM			0.30 Hrs	
198	Transfer	0.30	0.0	0.3 Hrs			4.54	0607:06	01:10 PM	0607:06	01:10 PM			01:10 PM	01:10 PM				
199	CIP	0.0	0.0	0.0 Hrs			4.55	0607:06	01:10 PM	0607:06	01:10 PM			01:10 PM	01:10 PM				
200	SIP	0.0	0.0	0.0 Hrs			4.55	0607:06	01:10 PM	0607:06	01:10 PM			01:10 PM	01:10 PM				
201	Clean Up	0.0	0.0	0.0 Hrs			4.55	0607:06	01:10 PM	0607:06	01:10 PM			01:10 PM	01:10 PM				
202	Subtotal	0.3			109.2		4.55	0607:06	01:10 PM	0607:06	01:10 PM								
203	9 B Homogenization																		
204	Set Up	0.00	0.0	0.0 Hrs	109.2		4.55	0607:06	01:10 PM	0607:06	01:10 PM	66.5 L @	1.6 LPM	01:51 PM	01:51 PM			0.68 Hrs	
205	Lysis	0.68	0.0	0.7 Hrs			4.55	0607:06	01:51 PM	0607:06	01:51 PM			01:51 PM	01:51 PM				
206	CIP	0.0	0.0	0.0 Hrs			4.58	0607:06	01:51 PM	0607:06	01:51 PM			01:51 PM	01:51 PM				
207	SIP	0.0	0.0	0.0 Hrs			4.58	0607:06	01:51 PM	0607:06	01:51 PM			01:51 PM	01:51 PM				
208	Clean Up	0.0	0.0	0.0 Hrs			4.58	0607:06	01:51 PM	0607:06	01:51 PM			01:51 PM	01:51 PM				
209	Subtotal	0.7			109.9		4.58	0607:06	01:51 PM	0607:06	01:51 PM								
210	10 B Heat Exchange																		
211	Set Up	0.50	0.0	0.5 Hrs	109.9		4.58	0607:06	01:51 PM	0607:06	01:51 PM	69.0 L @	3.8 LPM	02:09 PM	02:09 PM			0.30 Hrs	
212	Transfer	0.30	0.0	0.3 Hrs			4.58	0607:06	02:09 PM	0607:06	02:09 PM			02:09 PM	02:09 PM				
213	CIP	0.0	0.0	0.0 Hrs			4.59	0607:06	02:09 PM	0607:06	02:09 PM			02:09 PM	02:09 PM				
214	SIP	0.0	0.0	0.0 Hrs			4.59	0607:06	02:09 PM	0607:06	02:09 PM			02:09 PM	02:09 PM				
215	Clean Up	0.0	0.0	0.0 Hrs			4.59	0607:06	02:09 PM	0607:06	02:09 PM			02:09 PM	02:09 PM				
216	Subtotal	0.8			110.2		4.59	0607:06	02:09 PM	0607:06	02:09 PM								
217	8 C Heat Exchange																		
218	Set Up	0.00	0.0	0.0 Hrs	110.2		4.59	0607:06	02:09 PM	0607:06	02:09 PM	66.5 L @	3.7 LPM	02:27 PM	02:27 PM			0.30 Hrs	
219	Transfer	0.30	0.0	0.3 Hrs			4.59	0607:06	02:27 PM	0607:06	02:27 PM			02:27 PM	02:27 PM				
220	CIP	0.0	0.0	0.0 Hrs			4.60	0607:06	02:27 PM	0607:06	02:27 PM			02:27 PM	02:27 PM				
221	SIP	0.0	0.0	0.0 Hrs			4.64	0607:06	02:27 PM	0607:06	02:27 PM			02:27 PM	02:27 PM				
222	Clean Up	0.0	0.0	0.0 Hrs			4.69	0607:06	02:27 PM	0607:06	02:27 PM			02:27 PM	02:27 PM				
223	Subtotal	0.8			110.5		4.69	0607:06	02:27 PM	0607:06	02:27 PM								
224	9 C Homogenization																		
225	Set Up	0.00	0.0	0.0 Hrs	110.5		4.59	0607:06	02:09 PM	0607:06	02:09 PM	66.5 L @	3.7 LPM	02:27 PM	02:27 PM			0.30 Hrs	
226	Transfer	0.30	0.0	0.3 Hrs			4.59	0607:06	02:09 PM	0607:06	02:09 PM			02:27 PM	02:27 PM				
227	CIP	0.0	0.0	0.0 Hrs			4.60	0607:06	02:27 PM	0607:06	02:27 PM			02:27 PM	02:27 PM				
228	SIP	0.0	0.0	0.0 Hrs			4.64	0607:06	02:27 PM	0607:06	02:27 PM			02:27 PM	02:27 PM				
229	Clean Up	0.0	0.0	0.0 Hrs			4.69	0607:06	02:27 PM	0607:06	02:27 PM			02:27 PM	02:27 PM				
230	Subtotal	3.3			110.5		4.73	0607:06	02:27 PM	0607:06	02:27 PM								

FIG. 12.D

Process Time Line										Calculations	
Operation	Duration (Hrs.)		Rel. Time Scale (Hrs.)		Abs. Dvs.		Start		Finish		
	Calc.	A/D	Adj.	Prep	Exec.	Compl.	Start	End	Time	Date	
9 C Homogenization					13.5		06/07/96	06/07/96			
Set Up	0.00	0.0	0.0 Hrs				06/07/96	06/07/96	02:27 PM	06/07/96	
Lysis	0.60	0.0	0.7 Hrs	110.5	111.1		06/07/96	06/07/96	03:07 PM	06/07/96	66.5 L @ 1.6 LPM = 0.68 Hrs
CIP	1.0	0.0	1.0 Hrs			112.1	06/07/96	06/07/96	04:07 PM	06/07/96	
SIP	1.0	0.0	1.0 Hrs			113.1	06/07/96	06/07/96	05:07 PM	06/07/96	
Clean Up	1.0	0.0	1.0 Hrs			114.1	06/07/96	06/07/96	06:07 PM	06/07/96	
Sub Total	3.7		3.7 Hrs			114.1	06/07/96	06/07/96	06:07 PM	06/07/96	
10 C Heat Exchange											
Set Up	0.00	0.0	0.0 Hrs				06/07/96	06/07/96	03:07 PM	06/07/96	
Transfer	0.30	0.0	0.3 Hrs	111.1	111.4		06/07/96	06/07/96	03:25 PM	06/07/96	69.0 L @ 3.6 LPM = 0.30 Hrs
CIP	1.0	0.0	1.0 Hrs			112.4	06/07/96	06/07/96	04:25 PM	06/07/96	
SIP	1.0	0.0	1.0 Hrs			113.4	06/07/96	06/07/96	05:25 PM	06/07/96	
Clean Up	1.0	0.0	1.0 Hrs			114.4	06/07/96	06/07/96	06:25 PM	06/07/96	
Subtotal	3.3		3.3 Hrs			114.4	06/07/96	06/07/96	06:25 PM	06/07/96	
11 A Resolubilization											
Set Up	1.0	0.0	1.0 Hrs	108.9			06/07/96	06/07/96	11:52 AM	06/07/96	
Dilution	0.5	0.0	0.5 Hrs	109.4		4.49	06/07/96	06/07/96	12:52 PM	06/07/96	206.9 L @ 8.9 LPM = 0.50 Hrs
Agitate	0.5	0.0	0.5 Hrs	109.9		4.54	06/07/96	06/07/96	01:52 PM	06/07/96	0.50 Hrs
CIP	0.0	0.0	0.0 Hrs			4.56	06/07/96	06/07/96	01:52 PM	06/07/96	
SIP	0.0	0.0	0.0 Hrs			109.9	06/07/96	06/07/96	01:52 PM	06/07/96	
Clean Up	0.0	0.0	0.0 Hrs			109.9	06/07/96	06/07/96	01:52 PM	06/07/96	
Subtotal	2.0		2.0 Hrs			109.9	06/07/96	06/07/96	01:52 PM	06/07/96	
12 A Cent. Cent./Solids											
Set Up	1.0	0.0	1.0 Hrs	109.9			06/07/96	06/07/96	12:52 PM	06/07/96	
Centrifugation	0.5	0.0	0.5 Hrs	110.4		4.54	06/07/96	06/07/96	01:52 PM	06/07/96	275.9 L @ 9.2 LPM = 0.50 Hrs
Wash	0.1	0.0	0.1 Hrs	110.5		4.56	06/07/96	06/07/96	02:28 PM	06/07/96	3.0 L @ 0.5 LPM = 0.10 Hrs
CIP	0.0	0.0	0.0 Hrs			4.60	06/07/96	06/07/96	02:28 PM	06/07/96	8.0 L @ 0.5 LPM = 0.25 Hrs
SIP	0.0	0.0	0.0 Hrs			110.5	06/07/96	06/07/96	02:28 PM	06/07/96	
Clean Up	0.0	0.0	0.0 Hrs			110.5	06/07/96	06/07/96	02:28 PM	06/07/96	
Sub Total	1.6		1.6 Hrs			110.5	06/07/96	06/07/96	02:28 PM	06/07/96	
11 B Resolubilization											
Set Up	0.0	0.0	0.0 Hrs	110.5			06/07/96	06/07/96	02:28 PM	06/07/96	
Dilution	0.5	0.0	0.5 Hrs	111.0		4.59	06/07/96	06/07/96	03:13 PM	06/07/96	0.50 Hrs
Agitate	0.3	0.0	0.3 Hrs	111.2		4.62	06/07/96	06/07/96	03:49 PM	06/07/96	0.50 Hrs
CIP	1.0	0.0	1.0 Hrs			112.2	06/07/96	06/07/96	04:13 PM	06/07/96	
SIP	1.0	0.0	1.0 Hrs			113.2	06/07/96	06/07/96	05:13 PM	06/07/96	
Clean Up	1.0	0.0	1.0 Hrs			114.2	06/07/96	06/07/96	06:13 PM	06/07/96	
Subtotal	3.8		3.8 Hrs			114.2	06/07/96	06/07/96	06:13 PM	06/07/96	
12 B Cent. Cent./Solids											
Set Up	1.0	0.0	1.0 Hrs	111.2			06/07/96	06/07/96	02:13 PM	06/07/96	
Centrifugation	0.5	0.0	0.5 Hrs	111.7		4.59	06/07/96	06/07/96	03:43 PM	06/07/96	0.50 Hrs
Wash	0.1	0.0	0.1 Hrs	111.8		4.63	06/07/96	06/07/96	03:49 PM	06/07/96	0.10 Hrs
CIP	0.3	0.0	0.3 Hrs			4.66	06/07/96	06/07/96	04:04 PM	06/07/96	0.25 Hrs
SIP	1.0	0.0	1.0 Hrs			112.1	06/07/96	06/07/96	05:04 PM	06/07/96	
Clean Up	0.5	0.0	0.5 Hrs			113.1	06/07/96	06/07/96	05:34 PM	06/07/96	
Sub Total	3.4		3.4 Hrs			113.6	06/07/96	06/07/96	05:34 PM	06/07/96	
13 A Resolubilization											

FIG. 12E

Process Time Line										Rel. Time Scale (Hrs.)		Abs. Days		Start		Time		Calculations
Operation		Duration (Hrs.)		Calc. A/D		Adj.		Prep		Exec. Compl.		Start		Date		Time		
295										15.5			06/03/96	08:00 AM				
296	Set Up	1.0	0.0	1.0	0.0	110.5					4.56	06/07/96	01:28 PM	06/07/96	02:28 PM	60.7 L @ 2.0 LPM = 0.50 Hrs. 23.50 Hrs.		
297	Dilution	0.5	0.0	0.5	0.0	111.0					4.60	06/07/96	02:28 PM	06/07/96	02:58 PM			
298	Agitate	18.0	0.0	18.0	0.0	129.0					4.62	06/07/96	02:58 PM	06/03/96	08:58 AM			
299	CIP	1.0	0.0	1.0	0.0	130.0					5.37	06/07/96	08:58 AM	06/03/96	09:58 AM			
300	SIP	1.0	0.0	1.0	0.0	131.0					5.42	06/03/96	09:58 AM	06/03/96	10:58 AM			
301	Clean Up	1.0	0.0	1.0	0.0	132.0					5.46	06/03/96	10:58 AM	06/03/96	11:58 AM			
302	Subtotal	22.5		22.5		129.0					5.50	06/03/96	11:58 AM	06/03/96	01:58 PM			
303																		
304	14 A Concentration															28.99 SF		
305	Set Up	1.0	0.0	1.0	0.0	127.6					5.28	06/03/96	06:38 AM	06/03/96	07:38 AM	3.0 L/SF/Hr or 1.35 LPM 3.0 L/SF/Hr or 1.35 LPM 3.0 L/SF/Hr or 1.35 LPM 3.0 L/SF/Hr or 1.35 LPM 3.0 L/SF/Hr or 1.35 LPM		
306	Flush	0.7	0.0	0.7	0.0	128.3					5.32	06/03/96	07:38 AM	06/03/96	08:18 AM			
307	Prime	0.7	0.0	0.7	0.0	129.0					5.35	06/03/96	08:18 AM	06/03/96	08:58 AM			
308	Concentration	1.0	0.0	1.0	0.0	130.0					5.37	06/03/96	08:58 AM	06/03/96	09:58 AM			
309	Dilution	0.4	0.0	0.4	0.0	130.4					5.42	06/03/96	09:58 AM	06/03/96	10:25 AM			
310	Wash	0.9	0.0	0.9	0.0	131.3					5.43	06/03/96	10:25 AM	06/03/96	11:19 AM	3.0 L/SF/Hr or 1.35 LPM 3.0 L/SF/Hr or 1.35 LPM 3.0 L/SF/Hr or 1.35 LPM 3.0 L/SF/Hr or 1.35 LPM 3.0 L/SF/Hr or 1.35 LPM		
311	Flush	0.3	0.0	0.3	0.0	131.7					5.47	06/03/96	11:19 AM	06/03/96	12:19 PM			
312	Store	0.7	0.0	0.7	0.0	132.3					5.51	06/03/96	12:19 PM	06/03/96	01:19 PM			
313	CIP	1.0	0.0	1.0	0.0	133.3					5.56	06/03/96	01:19 PM	06/03/96	02:19 PM			
314	SIP	1.0	0.0	1.0	0.0	134.3					5.60	06/03/96	02:19 PM	06/03/96	03:19 PM			
315	Clean Up	1.0	0.0	1.0	0.0	135.3					5.64	06/03/96	03:19 PM	06/03/96	04:19 PM			
316	Sub Total	8.7		8.7		131.3					5.68	06/03/96	04:19 PM	06/03/96	05:19 PM	1.35 LPM		
317																		
318	15 A Microfiltration															12.60 SF		
319	Set Up	1.0	0.0	1.0	0.0	131.1					5.42	06/03/96	10:03 AM	06/03/96	11:03 AM	15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM		
320	Flush	0.1	0.0	0.1	0.0	131.2					5.46	06/03/96	11:03 AM	06/03/96	11:19 AM			
321	Prime	0.1	0.0	0.1	0.0	131.3					5.47	06/03/96	11:19 AM	06/03/96	11:49 AM			
322	Filtration	0.5	0.0	0.5	0.0	131.8					5.47	06/03/96	11:49 AM	06/03/96	11:51 AM			
323	Wash	0.0	0.0	0.0	0.0	131.8					5.49	06/03/96	11:51 AM	06/03/96	11:55 AM			
324	Regenerate	0.0	0.0	0.0	0.0	131.9					5.49	06/03/96	11:55 AM	06/03/96	12:55 PM	15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM		
325	Store	0.1	0.0	0.1	0.0	131.9					5.49	06/03/96	12:55 PM	06/03/96	01:55 PM			
326	CIP	1.0	0.0	1.0	0.0	132.9					5.54	06/03/96	01:55 PM	06/03/96	02:55 PM			
327	SIP	1.0	0.0	1.0	0.0	133.9					5.58	06/03/96	02:55 PM	06/03/96	03:55 PM			
328	Clean Up	1.0	0.0	1.0	0.0	134.9					5.62	06/03/96	03:55 PM	06/03/96	04:55 PM			
329	Sub Total	4.9		4.9		131.8					5.62	06/03/96	04:55 PM	06/03/96	05:55 PM	3.15 LPM		
330																		
331	16 A PIA MPLC															63.8 LCV		
332	Equilibration	1.1	0.0	1.1	0.0	131.4					5.43	06/03/96	10:17 AM	06/03/96	11:24 AM	100.0 CMHR or 4.76 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM		
333	Load	0.7	0.0	0.7	0.0	132.5					5.49	06/03/96	11:24 AM	06/03/96	12:31 PM			
334	Wash	1.3	0.0	1.3	0.0	133.9					5.52	06/03/96	12:31 PM	06/03/96	01:52 PM			
335	Elute A	0.0	0.0	0.0	0.0	135.2					5.58	06/03/96	01:52 PM	06/03/96	03:12 PM			
336	Elute B	0.0	0.0	0.0	0.0	135.2					5.63	06/03/96	03:12 PM	06/03/96	03:12 PM			
337	Regenerate	0.2	0.0	0.2	0.0	135.4					5.63	06/03/96	03:12 PM	06/03/96	03:25 PM	100.0 CMHR or 4.76 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM		
338	Store	0.4	0.0	0.4	0.0	135.9					5.64	06/03/96	03:25 PM	06/03/96	04:32 PM			
339	CIP	1.0	0.0	1.0	0.0	136.8					5.66	06/03/96	04:32 PM	06/03/96	05:32 PM			
340	SIP	1.0	0.0	1.0	0.0	137.9					5.70	06/03/96	05:32 PM	06/03/96	06:32 PM			
341	Clean Up	1.0	0.0	1.0	0.0	138.9					5.74	06/03/96	06:32 PM	06/03/96	07:32 PM			
342	Sub Total	5.2		5.2		135.2					5.79	06/03/96	07:32 PM	06/03/96	08:32 PM	4.76 LPM		
343																		
344	17 A PIA MPLC															Max FR		
345	Equilibration	0.6	0.0	0.6	0.0	135.6					5.62	06/03/96	02:59 PM	06/03/96	03:59 PM	4.76 LPM 100.0 CMHR or 1.58 LPM 50.0 CMHR or 0.79 LPM 50.0 CMHR or 0.79 LPM 50.0 CMHR or 0.79 LPM		
346	Load	1.1	0.0	1.1	0.0	136.3					5.63	06/03/96	03:59 PM	06/03/96	04:17 PM			
347	Wash	0.8	0.0	0.8	0.0	137.1					5.68	06/03/96	04:17 PM	06/03/96	05:03 PM			
348	Elute A	0.8	0.0	0.8	0.0	137.8					5.71	06/03/96	05:03 PM	06/03/96	05:49 PM			
349	Elute B	0.0	0.0	0.0	0.0	137.8					5.74	06/03/96	05:49 PM	06/03/96	06:52 PM			
350																34.75 CM Dia.		
351	Set Up	1.0	0.0	1.0	0.0	131.1					5.42	06/03/96	10:03 AM	06/03/96	11:03 AM	15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM		
352	Flush	0.1	0.0	0.1	0.0	131.2					5.46	06/03/96	11:03 AM	06/03/96	11:19 AM			
353	Prime	0.1	0.0	0.1	0.0	131.3					5.47	06/03/96	11:19 AM	06/03/96	11:49 AM			
354	Filtration	0.5	0.0	0.5	0.0	131.8					5.47	06/03/96	11:49 AM	06/03/96	11:51 AM			
355	Wash	0.0	0.0	0.0	0.0	131.8					5.49	06/03/96	11:51 AM	06/03/96	11:55 AM			
356	Regenerate	0.0	0.0	0.0	0.0	131.9					5.49	06/03/96	11:55 AM	06/03/96	12:55 PM	15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM 15.0 L/SF/Hr or 3.15 LPM		
357	Store	0.1	0.0	0.1	0.0	131.9					5.49	06/03/96	12:55 PM	06/03/96	01:55 PM			
358	CIP	1.0	0.0	1.0	0.0	132.9					5.54	06/03/96	01:55 PM	06/03/96	02:55 PM			
359	SIP	1.0	0.0	1.0	0.0	133.9					5.58	06/03/96	02:55 PM	06/03/96	03:55 PM			
360	Clean Up	1.0	0.0	1.0	0.0	134.9					5.62	06/03/96	03:55 PM	06/03/96	04:55 PM			
361	Sub Total	4.9		4.9		131.8					5.62	06/03/96	04:55 PM	06/03/96	05:55 PM	3.15 LPM		
362																		
363	18 A PIA MPLC															63.8 LCV		
364	Equilibration	1.1	0.0	1.1	0.0	131.4					5.43	06/03/96	10:17 AM	06/03/96	11:24 AM	100.0 CMHR or 4.76 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM		
365	Load	0.7	0.0	0.7	0.0	132.5					5.49	06/03/96	11:24 AM	06/03/96	12:31 PM			
366	Wash	1.3	0.0	1.3	0.0	133.9					5.52	06/03/96	12:31 PM	06/03/96	01:52 PM			
367	Elute A	0.0	0.0	0.0	0.0	135.2					5.58	06/03/96	01:52 PM	06/03/96	03:12 PM			
368	Elute B	0.0	0.0	0.0	0.0	135.2					5.63	06/03/96	03:12 PM	06/03/96	03:12 PM			
369	Regenerate	0.2	0.0	0.2	0.0	135.4					5.63	06/03/96	03:12 PM	06/03/96	03:25 PM	100.0 CMHR or 4.76 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM 50.0 CMHR or 2.38 LPM		
370	Store	0.4	0.0	0.4	0.0	135.9					5.64	06/03/96	03:25 PM	06/03/96	04:32 PM			
371	CIP	1.0	0.0	1.0	0.0	136.8					5.66	06/03/96	04:32 PM	06/03/96	05:32 PM			
372	SIP	1.0	0.0	1.0	0.0	137.9					5.70	06/03/96	05:32 PM	0				

FIG. 12F

Process Time Line		Rel. Time Scale (Hrs)				Act. Cntrs				Start				Finish				Calculations
		Calc.	Adj.	Prep	Exec.	Compl.	Start	End	Date	Time	Date	Time	Date	Time				
Operation																		
355	Regenerate	0.1	0.0	0.1 Hrs		135.0	5.74	5.75	06/03/98	05:00 AM	06/03/98	05:57 PM	06/03/98	05:57 PM	12.2 L @	100.0 CMHR	or	1.58 LPM
356	Store	0.3	0.0	0.3 Hrs		135.2	5.75	5.76	06/03/98	05:49 PM	06/03/98	06:13 PM	06/03/98	06:13 PM	24.4 L @	100.0 CMHR	or	1.58 LPM
357	CIP	1.0	0.0	1.0 Hrs		135.2	5.76	5.80	06/03/98	06:13 PM	06/03/98	07:13 PM	06/03/98	07:13 PM				
358	SIP	1.0	0.0	1.0 Hrs		140.2	5.80	5.84	06/03/98	07:13 PM	06/03/98	08:13 PM	06/03/98	08:13 PM				
359	Clean Up	1.0	0.0	1.0 Hrs		141.2	5.84	5.88	06/03/98	08:13 PM	06/03/98	09:13 PM	06/03/98	09:13 PM				
360	Sub Total	6.7		6.7 Hrs		137.8									Max FR			1.58 LPM
18 A Flow Dialysis															12.20 SF			
363	Set Up	1.0	0.0	1.0 Hrs	136.5		5.65	5.69	06/03/98	03:29 PM	06/03/98	04:29 PM	06/03/98	04:29 PM	24.4 L @	3.0 USFHR	or	0.81 LPM
364	Flush	0.7	0.0	0.7 Hrs	137.2		5.69	5.72	06/03/98	04:29 PM	06/03/98	05:09 PM	06/03/98	05:09 PM	24.4 L @	3.0 USFHR	or	0.81 LPM
365	Prime	0.7	0.0	0.7 Hrs	137.8		5.72	5.74	06/03/98	05:09 PM	06/03/98	05:49 PM	06/03/98	05:49 PM	36.0 L @	3.0 USFHR	or	0.81 LPM
366	Dialysis	1.0	0.0	1.0 Hrs	138.8		5.74	5.78	06/03/98	05:49 PM	06/03/98	06:49 PM	06/03/98	06:49 PM	0.0 L @	3.0 USFHR	or	0.81 LPM
367	Wash	0.0	0.0	0.0 Hrs	138.8		5.78	5.78	06/03/98	06:49 PM	06/03/98	07:09 PM	06/03/98	07:09 PM	12.2 L @	3.0 USFHR	or	0.81 LPM
368	Flush	0.3	0.0	0.3 Hrs		139.2	5.78	5.80	06/03/98	06:49 PM	06/03/98	07:09 PM	06/03/98	07:09 PM	24.4 L @	3.0 USFHR	or	0.81 LPM
369	Store	0.7	0.0	0.7 Hrs		139.8	5.80	5.83	06/03/98	07:09 PM	06/03/98	07:49 PM	06/03/98	07:49 PM				
370	CIP	1.0	0.0	1.0 Hrs		140.8	5.83	5.87	06/03/98	07:49 PM	06/03/98	08:49 PM	06/03/98	08:49 PM				
371	SIP	1.0	0.0	1.0 Hrs		141.8	5.87	5.91	06/03/98	08:49 PM	06/03/98	09:49 PM	06/03/98	09:49 PM				
372	Clean Up	1.0	0.0	1.0 Hrs		142.8	5.91	5.95	06/03/98	09:49 PM	06/03/98	10:49 PM	06/03/98	10:49 PM				
373	Sub Total	7.3		7.3 Hrs		138.8									Max FR			0.81 LPM
19 A PIA MPLC															7.0 L CV	0.4 HD		28.81 CM Dia.
376	Equilibration	0.5	0.0	0.5 Hrs	138.5		5.75	5.77	06/03/98	05:59 PM	06/03/98	06:31 PM	06/03/98	06:31 PM	34.8 L @	100.0 CMHR		1.09 LPM
377	Load	0.2	0.0	0.2 Hrs		139.1	5.78	5.79	06/03/98	06:49 PM	06/03/98	07:03 PM	06/03/98	07:03 PM	7.3 L @	50.0 CMHR		0.54 LPM
378	Wash	0.8	0.0	0.8 Hrs		139.7	5.79	5.82	06/03/98	07:40 PM	06/03/98	07:41 PM	06/03/98	07:41 PM	20.9 L @	50.0 CMHR		0.54 LPM
379	Elute A	0.6	0.0	0.6 Hrs		140.3	5.82	5.85	06/03/98	08:20 PM	06/03/98	08:20 PM	06/03/98	08:20 PM	20.9 L @	50.0 CMHR		0.54 LPM
380	Elute B	0.0	0.0	0.0 Hrs		140.3	5.85	5.85	06/03/98	08:20 PM	06/03/98	08:20 PM	06/03/98	08:20 PM	0.0 L @	30.0 CMHR		0.33 LPM
381	Regenerate	0.1	0.0	0.1 Hrs		140.4	5.85	5.86	06/03/98	08:20 PM	06/03/98	08:26 PM	06/03/98	08:26 PM	7.0 L @	100.0 CMHR		1.09 LPM
382	Store	0.2	0.0	0.2 Hrs		140.7	5.85	5.86	06/03/98	08:26 PM	06/03/98	08:39 PM	06/03/98	08:39 PM	13.9 L @	100.0 CMHR		1.09 LPM
383	CIP	1.0	0.0	1.0 Hrs		141.7	5.86	5.90	06/03/98	09:39 PM	06/03/98	10:39 PM	06/03/98	10:39 PM				
384	SIP	1.0	0.0	1.0 Hrs		142.7	5.90	5.94	06/03/98	10:39 PM	06/03/98	11:39 PM	06/03/98	11:39 PM				
385	Clean Up	1.0	0.0	1.0 Hrs		143.7	5.94	5.98	06/03/98	11:39 PM	06/03/98							
386	Sub Total	5.4		5.4 Hrs		140.3									Max FR			1.09 LPM
20 A Flow Dialysis															2.43 SF			
390	Set Up	0.0	0.0	0.0 Hrs	139.0		5.79	5.79	06/03/98	07:00 PM	06/03/98	07:00 PM	06/03/98	07:00 PM	4.9 L @	3.0 USFHR	or	0.12 LPM
391	Flush	0.7	0.0	0.7 Hrs	139.7		5.79	5.82	06/03/98	07:40 PM	06/03/98	08:20 PM	06/03/98	08:20 PM	4.9 L @	3.0 USFHR	or	0.12 LPM
392	Prime	0.7	0.0	0.7 Hrs	140.3		5.82	5.85	06/03/98	08:20 PM	06/03/98	10:20 PM	06/03/98	10:20 PM	14.8 L @	3.0 USFHR	or	0.12 LPM
393	Dialysis	2.0	0.0	2.0 Hrs		142.3	5.85	5.93	06/03/98	09:20 PM	06/03/98	10:20 PM	06/03/98	10:20 PM	0.0 L @	3.0 USFHR	or	0.12 LPM
394	Wash	0.0	0.0	0.0 Hrs	142.3		5.93	5.93	06/03/98	10:20 PM	06/03/98	10:40 PM	06/03/98	10:40 PM	2.4 L @	3.0 USFHR	or	0.12 LPM
395	Flush	0.3	0.0	0.3 Hrs		142.3	5.93	5.94	06/03/98	10:20 PM	06/03/98	11:20 PM	06/03/98	11:20 PM	4.9 L @	3.0 USFHR	or	0.12 LPM
396	Store	0.7	0.0	0.7 Hrs		143.3	5.94	5.97	06/03/98	11:20 PM	06/03/98	11:20 PM	06/03/98	11:20 PM				
397	CIP	0.0	0.0	0.0 Hrs		143.3	5.97	5.97	06/03/98	11:20 PM	06/03/98	11:20 PM	06/03/98	11:20 PM				
398	SIP	0.0	0.0	0.0 Hrs		143.3	5.97	5.97	06/03/98	11:20 PM	06/03/98	11:20 PM	06/03/98	11:20 PM				
399	Clean Up	0.0	0.0	0.0 Hrs		144.3	5.97	6.01	06/03/98	11:20 PM	06/03/98	12:20 AM	06/03/98	12:20 AM				
400	Sub Total	4.3		4.3 Hrs		142.3									Max FR			0.12 LPM
21 A PIA MPLC															6.3 L CV	0.4 HD		28.35 CM Dia.
404	Equilibration	0.5	0.0	0.5 Hrs	142.0		5.99	5.99	06/03/98	09:28 PM	06/03/98	09:57 PM	06/03/98	09:57 PM	26.8 L @	100.0 CMHR		0.81 LPM
405	Load	0.1	0.0	0.1 Hrs		142.4	5.99	5.99	06/03/98	10:28 PM	06/03/98	10:28 PM	06/03/98	10:28 PM	2.9 L @	50.0 CMHR		0.45 LPM
406	Wash	0.6	0.0	0.6 Hrs		143.0	5.99	5.99	06/03/98	11:31 PM	06/03/98	11:31 PM	06/03/98	11:31 PM	16.0 L @	50.0 CMHR		0.45 LPM
407	Elute A	0.6	0.0	0.6 Hrs		143.6	5.99	5.98	06/03/98	11:36 PM	06/03/98	11:36 PM	06/03/98	11:36 PM	0.0 L @	30.0 CMHR		0.27 LPM
408	Elute B	0.0	0.0	0.0 Hrs		143.6	5.98	5.98	06/03/98	11:36 PM	06/03/98	11:42 PM	06/03/98	11:42 PM	8.3 L @	100.0 CMHR		0.81 LPM
409	Regenerate	0.1	0.0	0.1 Hrs		143.7	5.98	5.99	06/03/98	11:36 PM	06/03/98	11:54 PM	06/03/98	11:54 PM	10.6 L @	100.0 CMHR		0.81 LPM
410	Store	0.2	0.0	0.2 Hrs		143.9	5.99	6.00	06/03/98	11:42 PM	06/03/98	11:54 PM	06/03/98	11:54 PM				
411	CIP	0.0	0.0	0.0 Hrs		143.9	6.00	6.00	06/03/98	11:54 PM	06/03/98							
412	SIP	0.0	0.0	0.0 Hrs		143.9	6.00	6.00	06/03/98	11:54 PM	06/03/98							
413	Clean Up	0.0	0.0	0.0 Hrs		143.9	6.00	6.00	06/03/98	11:54 PM	06/03/98							
414	Sub Total	0.0	0.0	0.0 Hrs		143.9	6.00	6.00	06/03/98	11:54 PM	06/03/98							

	Duration (Hrs.)		Rel. Time Scale (Hrs)			Abs. Days		Start		Finish		Calculations
	Calc.	AD	Asj.	Prep	Exec.	Comp.	Start	End	Date	Time	Time	
416	1.0	0.0	1.0 Hrs		15.5		06/03/98	06:00 AM	06/03/98	06:00 AM		Max FR 0.01 LPM
417	2.1		2.1 Hrs		143.8		06/03/98	11:54 PM	06/03/98	12:54 AM		
418												
419												2.2 L @ 50.0 USF/Hr or 0.07 LPM Max FR 0.07 LPM
420	0.5	0.0	0.5 Hrs	152.6			06/03/98	08:38 AM	06/03/98	08:38 AM		
421	0.5	0.0	0.5 Hrs	144.1			06/03/98	11:38 PM	06/03/98	12:38 AM		
422	0.5	0.0	0.5 Hrs		144.8		06/03/98	12:08 AM	06/03/98	12:38 AM		
423	0.0	0.0	0.0 Hrs		144.8		06/03/98	12:38 AM	06/03/98	12:38 AM		
424	0.0	0.0	0.0 Hrs		144.8		06/03/98	12:38 AM	06/03/98	12:38 AM		
425	1.0	0.0	1.0 Hrs		145.8		06/03/98	12:38 AM	06/03/98	01:38 AM		
426	1.5		1.5 Hrs		144.1							Max FR 0.07 LPM

FIG. 12H





FIG. 13




Solution Prep Vessel List/Procedure

Batch Tank		Batch Tank			Set Up Min.	Water Collect.		Weight/ Mix Min.	Ultrafiltration/Microfiltration				CIP			Total		Perc. Util.
No.	Min. LWV	No.	Min. LWV	Max. LWV		LPM	Min.		SF	L/SF/HR	Min.	Delay Factor	Adj. Min.	Cycle	Min.	SIP	Min.	
101	0.5	101	0.5	1	10		1	15	0.5	25	4.8	1.2	5.78			31.76	0.5	2%
102	1	102	1	2	10		1	15	1	25	4.8	1.2	5.78			31.76	0.5	
103	2	103	2	4	20		2	30	1	25	9.6	1.2	11.52			63.52	1.1	4%
104	4	104	4	10	20	10	1	30	2	25	12	1.2	14.4			65.4	1.1	
105	10	105	10	20	20	10	2	30	2	25	24	1.2	28.8			80.8	1.3	
106	20	106	20	50	20	10	5	30	10	25	12	1.2	14.4	CIP-1	60	109.4	1.8	3%
107	50	107	50	100	20	10	10	30	10	25	24	1.2	28.8	CIP-1	60	128.8	2.1	8%
108	100	108	100	250	0.5	50	5	30	30	25	20	1.2	24	CIP-1	60	99.5	1.7	18%
109	250	109	250	500	0.5	50	10	30	30	25	40	1.2	48	CIP-1	60	128.5	2.1	11%
110	500	110	500	1,500	1	50	30	30	60	25	60	1.2	72	CIP-1	60	173	2.9	10%
111	1500	111	1500	3,000	1	50	60	30	60	25	120	1.2	144	CIP-1	60	275	4.6	16%

1402 1406 1410 1404 1410 1412 1424 1506 1422 1434 1478

FIG. 15



Solution Prep Campaign Format

1626

Soln. ID	Storage Cond.		Soln. Prep Format			Solution Prep Cycles			Shelf Check	101 102 103 104 105 106 107						
	RT	4C	XP	MOD	BOD	BIA	Liters/ Batch	Preps/ Batch	Liters/ Prep	Days/ Bat. Cy.	Days/ Prep	Shelf Days	0.6	1	2	4
1 S-0101	X			0	X	0	1,666.50	1	1,666.50	7	7	56 OK				
2 S-0102	X			0	X	0	1.65	1	1.65	7	7	180 OK				
3 S-0103	X			0	X	0	1.65	1	1.65	7	7	180 OK				
4 S-0104	X			0	X	0	8.25	1	8.25	7	7	56 OK				
5 S-0105	X			0	X	0	8.25	1	8.25	7	7	56 OK				
6 S-0106	X			0	X	0	560.61	1	560.61	7	7	56 OK				
7 S-0107	X			0	X	0	125.93	1	125.93	7	7	56 OK				
8 S-0108	X			0	X	0	177.41	1	177.41	7	7	56 OK				
9 S-0109	X			0	X	0	22.18	1	22.18	7	7	56 OK				
10 S-0111	X			0	X	0	56.52	1	56.52	7	7	56 OK				
11 S-0112	X			0	X	0	113.03	1	113.03	7	7	56 OK				
12 S-0113	X			0	X	0	1,612.45	1	1,612.45	7	7	56 OK				
13 S-0114	X			0	X	0	574.10	1	574.10	7	7	56 OK				
14 S-0115	X			0	X	0	248.83	1	248.83	7	7	56 OK				
15 S-0116	X			0	X	0	497.65	1	497.65	7	7	56 OK				
16 S-0117	X			0	X	0	109.80	1	109.80	7	7	56 OK				
17 S-0118	X			0	X	0	497.65	1	497.65	7	7	56 OK				
18 S-0119	X			0	X	0	292.78	1	292.78	7	7	56 OK				
19 S-0120	X			0	X	0	109.80	1	109.80	7	7	56 OK				
20 S-0121	X			0	X	0	62.58	1	62.58	7	7	56 OK				
21 S-0122	X			0	X	0	0.00	1	0.00	7	7	56 OK				

1704 1610 1608 1620 1604 1614 1706

Fig. 17

Solution Prep Campaign Format

1626

Soln. ID	Tank Assignment				Solution Prep Schedule					Tank Flt						
	100 250	108 500	109 500	110 1500	111 4500 3000	Tank Assignment		Required By	Back Days	Avail. By	Hold Days	Inlt. Start	Float Days	Final Start	Next Prep	Prep. Hrs.
						Initial Assign.	Final Assign.									
1 S-0101					111	111	06/03/96	1	05/31/96	2	05/29/96	0	05/29/96	05/05/96	06/05/96	4
2 S-0102					111	102	08/05/96	1	08/04/96	2	05/31/96	0	05/31/96	05/31/96	06/07/96	4
3 S-0103						102	06/05/96	1	08/04/96	2	05/31/96	0	05/31/96	05/31/96	06/07/96	4
4 S-0104						104	06/05/96	1	06/04/96	2	05/31/96	0	05/31/96	05/31/96	06/07/96	4
5 S-0105						104	06/05/96	1	06/04/96	2	05/31/96	0	05/31/96	05/31/96	06/07/96	4
6 S-0106						110	06/07/96	1	08/06/98	2	08/04/98	0	06/04/98	06/04/98	06/11/98	4
7 S-0107				110		110	06/07/96	1	08/10/98	2	08/07/98	0	08/07/98	06/14/96	06/14/96	4
8 S-0108	108					108	08/11/98	1	08/10/98	2	08/07/98	0	08/07/98	06/14/96	06/14/96	4
9 S-0109	108					106	06/12/96	1	06/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
10 S-0111						107	06/12/96	1	06/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
11 S-0112	108					108	06/12/96	1	06/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
12 S-0113						111	06/12/96	1	08/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
13 S-0114				111		110	08/12/86	1	06/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
14 S-0115	108					108	08/12/86	1	06/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
15 S-0116						109	06/12/86	1	06/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
16 S-0117	108		109			108	06/12/96	1	06/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
17 S-0118			109			109	06/12/96	1	08/11/98	2	06/07/98	0	06/07/98	06/14/96	06/14/96	4
18 S-0119			109			109	06/12/96	1	06/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
19 S-0120	108					108	06/12/96	1	06/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
20 S-0121						107	06/12/96	1	06/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4
21 S-0122						0	06/12/96	1	08/11/96	2	06/07/96	0	06/07/96	06/14/96	06/14/96	4

Min. 05/29/96 08/14/96

Min. 08/03/96 06/12/96


Min Max Sat Sun

1722 1726 1728

1724

FIG 18

00000000000000000000



306




FIG. 20

2102

2104

2106

	Category/Assay	Code	Man Hour			Disp. Material
			Set Up	Per Sample	Clean Up	
1	Environmental					
2	Temperature	E-1	0.5	0.1	0.5	
3	Humidity	E-2	0.5	0.1	0.5	
4	Particle Count	E-3	0.5	0.2	0.5	
5						
6	Analytical					
7	Visual					
8	Certificate of Analysis	AV-1	0.25	0.2	0.5	
9	Appearance	AV-2	0.25	0.05	0.25	
10	Chemical					
11	Solubility	AC-1	0.5	0.1	0.5	
12	pH	AC-2	0.25	0.05	0.25	
13	Osmolality	AC-3	0.25	0.1	0.25	
14	Water Content (by Karl Fischer)	AC-4	0.5	0.2	0.5	
15	Key Element Analysis (by ICP Atomic Adsorption Spectroscopy)	AC-5	1	0.25	1	
16	GC/Mass Spec	AC-6	1	0.25	1	
17	Biochemical					
18	DNA					
19	DNA Fluorochrome Stain	AB-1	0.5	0.1	0.5	
20	Protein					
21	Hemoglobin	AB-2	0.5	0.1	0.5	
22	Electrophoretic Profiles by SDS-PAGE	AB-3	1	0.2	1	
23	A280	AB-4	0.25	0.1	0.25	
24	Bradford Assay	AB-5	0.5	0.1	0.5	
25	Amino Acid Analysis by HPLC	AB-6	1	0.25	1	
26	Endotoxin		0.5	0.1	0.5	
27	Gel Clot LAL	AB-7				
28	Immunological					
29	ELISA	AI-1	1	0.1	1	
30	Western Blots	AI-2	1.5	0.2	1.5	
31	Activity					
32	Chromagenic Substrate Assays	AA-1	1	0.1	1	
33						
34	In Vitro Biological					
35	Microbiological	VB-1	0.5	0.2	0.5	
36	Mycoplasma (Barile Method)	VB-2	0.5	0.2	0.5	
37	Bacteriophage (Screened)	VB-3	0.5	0.2	0.5	
38	Cell Passage Test	VB-4	1	0.2	1	
39	Adventitious viral Agents		2	0.2	1	
40	CPE	VB-5	2	0.2	1	
41	BVD	VB-6	2	0.2	1	
42	P13	VB-7	2	0.2	1	
43	IBR	VB-8	2	0.2	1	
44	Virus Neutralization Titers (9CFR)					
45	BVD	VB-9	2	0.2	1	
46	P13	VB-10	2	0.2	1	
47	IBR	VB-11	2	0.2	1	
48	Tritiated Thymidine Uptake in Mouse Cells	VB-12	2	0.2	1	
49	General Safety Test (Guinea Pigs)	VB-13	1	0.2	1	
50						
51						

FIG. 21




FIG. 22






FIG. 23






FIG. 26

2025 RELEASE UNDER E.O. 14176




FIG. 27




FIG. 28












FIG. 29

1316

1318




FIG. 30

2804

Q110 7906Z

2102

31

1320





FIG. 32

FIG. 33



3302




FIG. 34

3304




FIG. 35

Prep Equipment Protocol - Bench Sink

3408

3602

3604

	Cycler Code	Minutes/Cycle									Total	
		Load	Pre Wash Rinse		Detergent Wash			Post Wash Rinse		Final Rinse		Hold/ Dry
			NPHW	NPCW	Minutes	Reagent	Gm/CF	NPHW	NPCW			
1	B8-1	5	2	2	5	Alconox	0.5	2	2	2		20
2	B8-2	5	2	2	5	Alconox	0.5	2	2	2		20
3	B8-3	5	2	2	5	Alconox	0.5	2	2	2		20
4	B8-4	5	2	2	5	Alconox	0.5	2	2	2		20
5	B8-5	5	2	2	5	Alconox	0.5	2	2	2		20

FIG. 36A

3408

Protocol		Minutes/Cycle								Total
Cycle Code	Load	Pre Wash Rinse		Detergent Wash			Post Wash Rinse		Final Rinse	
		NPHW	NPCW	Minutes	Reagent	Gm/CF	NPHW	NPCW		
1	WS-1	5	2	2	5 Alconox	0.5	2	2	2	
2	WS-2	5	2	2	6 Alconox	0.5	2	2	2	
3	WS-3	5	2	2	5 Alconox	0.5	2	2	2	
4	WS-4	5	2	2	5 Alconox	0.5	2	2	2	
5	WS-5	5	2	2	5 Alconox	0.5	2	2	2	

1. The first part of the document is a list of references. The references are listed in a standard format, with the author's name, the title of the work, and the publisher. The references are as follows:

1. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
2. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
3. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
4. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
5. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
6. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
7. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
8. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
9. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
10. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.

2. The second part of the document is a list of references. The references are listed in a standard format, with the author's name, the title of the work, and the publisher. The references are as follows:

1. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
2. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
3. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
4. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
5. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
6. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
7. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
8. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
9. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.
10. J. H. Van Veen, "Acoustic signal processing for hearing aids," *IEEE Transactions on Speech and Audio Processing*, vol. 1, no. 1, pp. 3-14, 1993.

FIG. 36B

3408

	Cycle Code	Minutes/Cycle										Total
		Load	Pre Wash Rinse		Detergent Wash			Post Wash Rinse		Final Rinse	Unload	
			NPHW	NPCW	Minutes	Reagent	Gm/CF	NPHW	NPCW			
1	GW-1	15	2	2	5	Alconox	0.5	2	2	2	10	40
2	GW-2	15	2	2	6	Alconox	0.5	2	2	2	10	40
3	GW-3	15	2	2	6	Alconox	0.5	2	2	2	10	40
4	GW-4	15	2	2	5	Alconox	0.5	2	2	2	10	40
5	GW-5	15	2	2	5	Alconox	0.5	2	2	2	10	40

FIG. 36C

Prep Equipment Protocol - Glassware Dryer

3408

	Cycle Code	Load	Heat Up Minutes	Dry		Cool Minutes	Unload	Total
				Temp (C)	Minutes			
1	DO-1	10	30	250	40	30	10	120
2	DO-2	10	30	250	25	30	10	105
3	DO-3	10	30	250	25	30	10	105
4	DO-4	10	30	250	25	30	10	105
5	DO-5	10	30	250	25	30	10	105

3618

3620

3622

3624

3626

3628

FIG. 36D

666730" 00000000

Prep Equipment Protocol - Carboy Washer

3408

1.55 min

sum

Minutes/Cycle										
Load	Pre Wash Rinse		Detergent			Post Wash Rinse		Final Rinse	Unload	Total
	NPHW	NPCW	Minutes	Reagent	Gm/CF	NPHW	NPCW			
15	2	2	5	Alconox	0.5	2	2	2	15	15
15	2	2	5	Alconox	0.5	2	2	2	15	15
15	2	2	5	Alconox	0.5	2	2	2	15	15
15	2	2	5	Alconox	0.5	2	2	2	15	15
15	2	2	5	Alconox	0.5	2	2	2	15	15
15	2	2	5	Alconox	0.5	2	2	2	15	15

FIG. 36E

3408

	Cycle Code	Load	Heat Up Minutes	Dry		Cool Minutes	Unload	Total
				Temp (C)	Minutes			
1	CD-1	10	30	250	40	30	10	100
2	CD-2	10	30	250	25	30	10	85
3	CD-3	10	30	260	25	30	10	85
4	CD-4	10	30	260	25	30	10	85
5	CD-5	10	30	250	25	30	10	85

FIG. 36F

1. The first part of the document discusses the importance of maintaining accurate records of all transactions and activities related to the business. It emphasizes the need for transparency and accountability in financial reporting.

2. The second part of the document outlines the various methods and techniques used to collect and analyze data. It includes a detailed description of the experimental setup and the procedures followed to ensure the reliability of the results.

3. The third part of the document presents the results of the study, showing the trends and patterns observed in the data. It includes a comparison of the findings with previous research and a discussion of the implications for future work.

4. The fourth part of the document concludes the study, summarizing the key findings and providing recommendations for further research. It also includes a list of references and a list of figures and tables.

3606 3608 3610 3612v Prep Equipment Protocol - Steam Sterilizer

3614 3616

3606 3608

Cycles										SS-1					SS-2					SS-3				
	Press. (Bar)	Minutes To Ach.	Minutes To Hold	No. of Cycles	Subt.	Press. (Bar)	Minutes To Ach.	Minutes To Hold	No. of Cycles	Subt.	Press. (Bar)	Minutes To Ach.	Minutes To Hold	No. of Cycles	Subt.	Press. (Bar)	Minutes To Ach.	Minutes To Hold	No. of Cycles	Subt.				
1										20											20			
2																								
3																								
4			16	1	16																			
5																								
6																								
7																								
8					18																			
9																								
10																								
11																								
12																								
13					60																			
14																								
15																								
16																								
17																								
18																								
19																								
20																								
21																								
22																								
23																								
24																								
25																								
26																								
27																								
28																								
29																								
30																								
31																								
32					161					196											230			
33					2.7					3.3											3.8			

FIG. 365

Prep Equipment Protocol - Dry Heat Sterilizer

340°

	Cycle Code	Load	Heat Up Minutes	Sterilization		Cool Minutes	Unload	Total
				Temp (C)	Minutes			
1	SO-1	15	30	250	40	30	15	130
2	SO-2	15	30	260	25	30	15	115
3	SO-3	15	30	250	25	30	15	115
4	SO-4	15	30	250	25	30	15	115
5	SO-5	15	30	260	25	30	15	115

FIG. 36 H

20250303 09:00:00

Prep-Equipment Protocol - Equipment Prep Procedures

3706

3702

3704

			EPC1	EPC2	EPC3	EPC4	EPC5	EPC6	EPC7
1	Initial Rinse								
2									
3	Bench Sink - 1								
4	Procedure Protocol		BS-1	BS-1	BS-2	BS-1			
5	Duration	PHrs.	0.33	0.33	0.33	0.33			
6	Hold/Dry	PHrs.	0	0	0				
7	Subtotal	PHrs.	0.33	0.33	0.33	0.33	0.00	0.00	0.00
8	Cummulative	PHrs.	0.33	0.33	0.33	0.33	0.00	0.00	0.00
9									
10	Wash Station - 1								
11	Procedure Protocol						WS-1	WS-1	
12	Duration	PHrs.					0.25	0.25	
13	Hold/Dry	PHrs.							
14	Subtotal	PHrs.	0.00	0.00	0.00	0.00	0.25	0.25	0.00
15	Cummulative	PHrs.	0.33333	0.33333	0.33333	0.33333	0	0	0
16									
17	Cleaning								
18									
19	Bench Sink - 1								
20	Procedure Protocol		BS-3	BS-3	BS-4				
21	Duration	PHrs.	0.33	0.33	0.33				
22	Hold/Dry	PHrs.							
23	Subtotal	PHrs.	0.33	0.33	0.33	0.00	0.00	0.00	0.00
24	Cummulative	PHrs.	0.66667	0.66667	0.66667	0.33333	0	0	0
25									
26	Glassware Washer - 1								
27	Procedure Protocol					GW-1			
28	Duration	PHrs.				0.67			
29	Hold/Dry	PHrs.							
30	Subtotal	PHrs.	0.00	0.00	0.00	0.67	0.00	0.00	0.00
31	Cummulative	PHrs.	0.66667	0.66667	0.66667	1	0	0	0
32									
33	Glassware Dryer - 1								
34	Procedure Protocol		GD-1	GD-1	GD-2	GD-3			
35	Duration	PHrs.	2.00	2.00	1.75	1.75			
36	Hold/Dry	PHrs.							
37	Subtotal	PHrs.	2.00	2.00	1.75	1.75	0.00	0.00	0.00
38	Cummulative	PHrs.	2.66667	2.66667	2.41667	2.75	0	0	0
39									
40	Carboy Washer - 1								
41	Procedure Protocol						CW-1	CW-1	
42	Duration	PHrs.					0.25	0.25	
43	Hold/Dry	PHrs.							
44	Subtotal	PHrs.	0.00	0.00	0.00	0.00	0.25	0.25	0.00
45	Cummulative	PHrs.	2.66667	2.66667	2.41667	2.75	0.25	0.25	0
46									
47	Carboy Dryer - 1								
48	Procedure Protocol						CD-1	CD-1	
49	Duration	PHrs.					1.67	1.67	
50	Hold/Dry	PHrs.							
51	Subtotal	PHrs.	0.00	0.00	0.00	0.00	1.67	1.67	0.00
52	Cummulative	PHrs.	2.66667	2.66667	2.41667	2.75	1.91667	1.91667	0
53									
54	Prep								
55									
56	Staffing		2	2	2	2	2	2	2
57									
58	Preassembly								
59	Man Hours	MHrs.		1					
60	Procedure Hours			0.5					

FIG. 37A

Prep Equipment Protocol - Equipment Prep Procedures

			EPC1	EPC2	EPC3	EPC4	EPC5	EPC6	EPC7
61	Cummulative	PHrs.	2.68667	3.16667	2.41667	2.75	1.91667	1.91667	0
62									
63	Wrap								
64	Man Hours	MHrs.	1.5	1.5	1.5	1.5	1.5	1.5	1.5
65	Procedure Hours		0.75	0.75	0.75	0.75	0.75	0.75	0.75
66	Cummulative	PHrs.	3.41667	3.91667	3.16667	3.5	2.66667	2.66667	0.75
67									
68	Sterilization								
69									
70	Autoclave - 1								
71	Procedure		SS-1	SS-1	SS-1	SS-1	SS-2		SS-3
72	Duration	PHrs.	2.68	2.68	2.68	2.68	3.25		3.83
73	Hold/Dry	PHrs.							
74	Subtotal	PHrs.	2.68	2.68	2.68	2.68	3.25	0.00	3.83
75	Cummulative	PHrs.	6.10	6.60	5.85	6.18	5.92	2.67	4.58
76									
77	Dry Heat - 1								
78	Procedure							SO-1	
79	Hours/Load	PHrs.						2.17	
80	Hold/Dry	PHrs.							
81	Subtotal	PHrs.	0.00	0.00	0.00	0.00	0.00	2.17	0.00
82	Cummulative	PHrs.	6.10	6.60	5.85	6.18	5.92	4.83	4.58
83									
84	Total		6.10	6.60	5.85	6.18	6.17	5.08	4.58
85									
86	Max		2.68	2.68	2.68	2.68	3.25	2.17	3.83

FIG. 378

3306




FIG. 38

3902

3904 - 3906
3908
3910 3912 3914 3916 3918 3920 3922

53
FH

3308




FIG. 40

3310




FIG. 41

Equipment Prep Load Table

TASK -Equipment Name-	Unit Oper End Time		EPC-1		EPC-2										EPC-3	
	Date	Time	Specialty Glass Siphon Tubes	Total	Instruments-			Fittings			Reducers	Hose Barbs	Clamps	Total CF	Plasticware Beakers	
					PI 0.03	DO Probe 0.08	pH Probe 0.06	Tees 0.03	Elbows 0.02	Crosses 0.06						
1 Inoculum Prep	06/04/98	02:30 PM		0										0.00		
2 Flask Growth	06/05/98	01:30 PM		0										0.00		
3 Seed Fermentation	06/06/98	03:30 PM		0										0.00		
4 Fermentation	06/07/98	12:00 PM		0	4 0.111			6 0.17			2 0.03	4 0.03	16 0.17	0.50		
5 Heat Exchange	06/07/98	01:00 PM		0	3			4 0.11				4 0.03	8 0.08	0.22		
6 Cont. Cent/Solids	06/07/98	11:51 AM		0	3 0.083			4 0.11				4 0.03	8 0.08	0.31		
1 Inoculum Prep	06/08/98	02:30 PM		0										0.00		
2 Flask Growth	06/07/98	01:30 PM		0										0.00		
3 Seed Fermentation	06/08/98	03:30 PM		0										0.00		
4 Fermentation	06/09/98	09:00 AM		0	4 0.111			6 0.17			2 0.03	4 0.03	16 0.17	0.50		
5 Heat Exchange	06/09/98	10:00 AM		0	3 0.083			4 0.11				4 0.03	8 0.08	0.31		
6 Cont. Cent/Solids	06/09/98	08:51 AM		0	3 0.083			4 0.11				4 0.03	8 0.08	0.31		
1 Inoculum Prep	06/08/98	02:30 PM		0										0.00		
2 Flask Growth	06/09/98	01:30 PM		0										0.00		
3 Seed Fermentation	06/10/98	03:30 PM		0										0.00		
4 Fermentation	06/03/98	10:00 AM		0	4 0.111			6 0.17			2 0.03	4 0.03	16 0.17	0.50		
5 Heat Exchange	06/11/98	09:00 AM		0	3 0.083			4 0.11				4 0.03	8 0.08	0.31		
6 Cont. Cent/Solids	06/11/98	08:51 AM		0	3 0.083			4 0.11				4 0.03	8 0.08	0.31		
1 Inoculum Prep	06/11/98	12:15 PM		0										0.00		
2 Cell Resuspension	06/11/98	09:33 AM		0										0.00		
3 Heat Exchange	06/11/98	09:51 AM		0										0.00		
4 Cell Disruption	06/11/98	10:09 AM		0										0.00		
5 Heat Exchange	06/11/98	10:09 AM		0										0.00		

4210

[illegible]

FFs: 42B

06/11/96 06/12/96 06/13/96

Equipment Prep Load Table

4218

4220

Equipment Items	Unit Oper End Time		EPC-1		EPC-2										EPC-3	
			Specialty Glass		Total	Instruments				Fittings			Plasticware			
	Date	Time	Siphon Tubes	PI 0.03		DO Probe 0.06	pH Probe 0.06	Tees 0.03	Elbows 0.02	Crosses 0.08	Reducers 0.01	Hose Barbs 0.01	Clamps 0.01	Total CF	Plasticware Beakers 0.03	
8 Heat Exchange	06/11/96	10:27 AM			0										0.00	
9 Cell Disruption	06/11/96	10:45 AM			0										0.00	
10 Heat Exchange	06/11/96	12:00 AM			0										0.00	
8 Heat Exchange	06/11/96	02:21 PM			0										0.00	
9 Cell Disruption	06/11/96	02:39 PM			0										0.00	
10 Heat Exchange	06/11/96	02:57 PM			0										0.00	
11 IB Resuspension	06/11/96	10:57 AM			0										0.00	
12 Centrifugation	06/11/96	11:33 AM			0										0.00	
11 IB Resuspension	06/11/96	03:08 PM			0										0.00	
12 Centrifugation	06/11/96	03:12 PM			0										0.00	
13 Renaturation	06/12/96	08:43 AM			0										0.00	
14 Buffer Exchange	06/12/96	11:47 AM			0										0.00	
15 Clarification	06/12/96	11:03 AM			0										0.00	
16 Chromatography 1	06/12/96	03:59 PM			0										0.00	
17 Chromatography 2	06/12/96	06:59 PM			0										0.00	
18 Buffer Exchange	06/12/96	08:27 PM			0										0.00	
19 Chromatography 3	06/12/96	10:07 PM			0										0.00	
20 Buffer Exchange	06/12/96	10:38 PM			0										0.00	
21 Chromatography 4	06/13/96	12:14 AM			0										0.00	
22 Sterile Filtration	06/13/96	12:48 AM			0										0.00	
Totals															3.25	

FIG. 42C

4222

4224


4226

4228

Equipment, Top Load Table

Equipment Items	Unit Oper End Time DateTime		EPC-4										EPC-5						EPC-6					
			Rubber Stoppers		Flexible Tubing		Total CF	Small Glassware		Total CF	PP Carboys			Total CF	BSG Carboys			Total CF						
			Flasks 0.25	Silicone 0.00	Butyl 0.03	Silicone 0.33		Neoprene 3.33	Beakers 0.03		Flasks 0.25	10L 1.3333	20L 4.88		45L 10.7	10L 1.3333	20L 4.88		45L 10.7					
8 Heat Exchange	08/11/96	10:27 AM						0.00			0									0.00				
9 Cell Disruption	08/11/96	10:45 AM						0.00			0									0.00				
10 Heat Exchange	08/11/96	12:00 AM						0.00			5 1.25	1.25								0.00				
8 Heat Exchange	08/11/96	02:21 PM						0.00				0								0.00				
9 Cell Disruption	06/11/98	02:39 PM						0.00				0								0.00				
10 Heat Exchange	08/11/96	02:57 PM						0.00			5 1.25	1.25								0.00				
11 IB Resuspension	08/11/96	10:57 AM						0.00				0								0.00				
12 Centrifugation	08/11/96	11:33 AM						0.00				0								0.00				
11 IB Resuspension	08/11/96	03:06 PM						0.00				0								0.00				
12 Centrifugation	08/11/96	03:12 PM						0.00				0								0.00				
13 Renaturation	08/12/98	08:43 AM						0.00				0								0.00				
14 Buffer Exchange	08/12/96	11:47 AM						0.00				0								0.00				
15 Clarification	08/12/96	11:03 AM						0.00				0								0.00				
16 Chromatography 1	08/12/96	03:59 PM						0.00				0								0.00				
17 Chromatography 2	08/12/96	06:59 PM						0.00				0								0.00				
18 Buffer Exchange	08/12/98	08:27 PM						0.00				0								0.00				
19 Chromatography 3	08/12/96	10:07 PM						0.00				0								0.00				
20 Buffer Exchange	08/12/98	10:38 PM						0.00				0								0.00				
21 Chromatography 4	08/13/96	12:14 AM						0.00				0								0.00				
22 Sterile Filtration	08/13/98	12:48 AM						0.00				0								0.00				
Totals																								

FIG. 42 D






FIG. 44

4506

QC Load Table - PE Module

4504

4502

	Operation	QA/QC Samples										Immunochemical	Act.												
		Start		Finish		Visual	Chemical							Biological											
		Date	Time	Date	Time		AV-1	AV-2	AC-1	AC-2	AC-3			AC-4	AC-5	AC-6	AB-1	AB-2	AB-3	AB-4	AB-5	AB-6	AB-7	AI-1	AI-2
1	1 A Inoculum Prep	06/03/96	08:00 AM																						
2																									
3	Set Up	06/03/96	09:30 AM	06/03/96	12:30 PM																				
4	Preincubation	06/03/96	12:30 PM	06/03/96	03:30 PM																				
5	Incubation	06/03/96	03:30 PM	06/04/96	02:30 PM																				
6	Clean Up	06/04/96	02:30 PM	06/04/96	02:45 PM																				
6	Subtotal																								
7																									
8	2 A Flask Growth																								
9																									
10	Set Up	06/04/96	12:30 PM	06/04/96	01:30 PM																				
11	Preincubation	06/04/96	01:30 PM	06/04/96	02:30 PM																				
12	Incubation	06/04/96	02:30 PM	06/05/96	01:30 PM																				
13	Clean Up	06/05/96	01:30 PM	06/05/96	01:45 PM																				
13	Subtotal																								
14																									
16	3 A Seed Fermentation																								
17	Set Up	06/05/96	11:30 AM	06/05/96	12:30 PM																				
18	Preincubation	06/05/96	12:30 PM	06/05/96	01:30 PM																				
19	Fermentation	06/05/96	01:30 PM	06/06/96	10:30 AM																				
20	Harvest	06/06/96	10:30 AM	06/06/96	11:00 AM																				
21	CIP	06/06/96	10:30 AM	06/06/96	11:30 AM																				
22	SIP	06/06/96	11:30 AM	06/06/96	12:30 PM																				
23	Clean Up	06/06/96	12:30 PM	06/06/96	03:30 PM																				
24	Subtotal																								
25																									
26	4 A Production Fermentation																								
27	Set Up	06/06/96	09:00 AM	06/06/96	10:00 AM																				
28	Preincubation	06/06/96	10:00 AM	06/06/96	11:00 AM																				
29	Fermentation	06/06/96	11:00 AM	06/07/96	08:00 AM																				
30	CIP	06/07/96	08:00 AM	06/07/96	09:00 AM																				
31	SIP	06/07/96	09:00 AM	06/07/96	10:00 AM																				
32	Clean Up	06/07/96	10:00 AM	06/07/96	12:00 PM																				
33	Subtotal																								
34																									
35																									
36	6 A Heat Exchange																								
37	Set Up	06/07/96	08:00 AM	06/07/96	08:30 AM																				
38	Transfer	06/07/96	08:30 AM	06/07/96	09:00 AM																				
39	CIP	06/07/96	09:00 AM	06/07/96	10:00 AM																				
40	SIP	06/07/96	10:00 AM	06/07/96	11:00 AM																				
41	Clean Up	06/07/96	11:00 AM	06/07/96	01:00 PM																				
42	Subtotal																								
43																									
44																									
45	6 A Cont. Cent./Solids																								
46	Set Up	06/07/96	08:00 AM	06/07/96	09:00 AM																				
47																									

FIG. 45A

4506

QC Load Table - PE Module

Operation		Q/JOC Samples										Biochemical						Immunological		
		Start		Finish		Visual		Chemical						Biochemical						Act.
		Date	Time	Date	Time	AV-1	AV-2	AC-1	AC-2	AC-3	AC-4	AC-5	AC-6	AB-1	AB-2	AB-3	AB-4	AB-5	AB-6	
48	Centrifugation	06/03/98	08:00 AM	06/07/98	10:00 AM											3				3
49	Wash	06/07/98	09:00 AM	06/07/98	10:06 AM															
50	CIP	06/07/98	10:06 AM	06/07/98	10:21 AM															
51	SIP	06/07/98	10:21 AM	06/07/98	11:21 AM															
52	Clean Up	06/07/98	11:21 AM	06/07/98	11:51 AM															
53	Sub Total																			
54																				
55	1 B Inoculum Prep																			
56	Set Up	06/03/98	01:30 PM	06/03/98	02:30 PM															
57	Preincubation	06/03/98	02:30 PM	06/03/98	03:30 PM															
58	Incubation	06/03/98	03:30 PM	06/04/98	02:30 PM															
59	Clean Up	06/04/98	02:30 PM	06/04/98	02:45 PM															
60	Subtotal																			
61																				
62	2 B Flask Growth																			
63	Set Up	06/04/98	12:30 PM	06/04/98	01:30 PM															
64	Preincubation	06/04/98	01:30 PM	06/04/98	02:30 PM															
65	Incubation	06/04/98	02:30 PM	06/05/98	01:30 PM															
66	Clean Up	06/05/98	01:30 PM	06/05/98	01:45 PM															
67	Subtotal																			
68																				
69	3 B Seed Fermentation																			
70	Set Up	06/05/98	11:30 AM	06/05/98	12:30 PM															
71	Preincubation	06/05/98	12:30 PM	06/05/98	01:30 PM															
72	Fermentation	06/05/98	01:30 PM	06/06/98	10:30 AM															
73	Harvest	06/06/98	10:30 AM	06/06/98	11:00 AM															
74	CIP	06/06/98	10:30 AM	06/06/98	11:30 AM															
75	SIP	06/06/98	11:30 AM	06/06/98	12:30 PM															
76	Clean Up	06/06/98	12:30 PM	06/06/98	03:30 PM															
77	Subtotal																			
78																				
79	4 B Production Fermentation																			
80	Set Up	06/06/98	09:00 AM	06/06/98	10:00 AM															
81	Preincubation	06/06/98	10:00 AM	06/06/98	11:00 AM															
82	Fermentation	06/06/98	11:00 AM	06/07/98	08:00 AM															
83	CIP	06/07/98	08:00 AM	06/07/98	08:00 AM															
84	SIP	06/07/98	09:00 AM	06/07/98	10:00 AM															
85	Clean Up	06/07/98	10:00 AM	06/07/98	12:00 PM															
86	Subtotal																			
87																				
88	6 B Heat Exchange																			
89	Set Up	06/07/98	08:00 AM	06/07/98	08:30 AM															
90	Transfer	06/07/98	08:00 AM	06/07/98	09:00 AM															
91	CIP	06/07/98	08:00 AM	06/07/98	10:00 AM															
92																				
93																				
94																				

FIG- 45B

[illegible]

Fig. 45c

Fig. 45c

Fig. 45D

[illegible]

Fig. 45E

QC Load Table - PE Module

06/07/96 06:00 AM

	Operation	QA/QC Samples										Immunological	Act.										
		Start		Finish		Chemical								Biochemical									
		Date	Time	Date	Time	AV-1	AV-2	AC-1	AC-2	AC-3	AC-4			AC-5	AC-6	AB-1	AB-2	AB-3	AB-4	AB-5	AB-6	AB-7	AI-1
240	CIP	06/03/98	08:00 AM	06/07/98	04:07 PM																		
241	SIP	06/07/98	04:07 PM	06/07/98	05:07 PM																		
242	Clean Up	06/07/98	05:07 PM	06/07/98	06:07 PM																		
243	Sub Total																						
244																							
245	10 C Heat Exchange																						
246	Set Up	06/07/98	03:07 PM	06/07/98	03:07 PM																		
247	Transfer	06/07/98	03:07 PM	06/07/98	03:25 PM																		
248	CIP	06/07/98	03:25 PM	06/07/98	04:25 PM																		
249	SIP	06/07/98	04:25 PM	06/07/98	05:25 PM																		
250	Clean Up	06/07/98	05:25 PM	06/07/98	06:25 PM																		
251	Subtotal																						
252																							
253																							
254	11 A Resolubilization																						
255	Set Up	06/07/98	11:52 AM	06/07/98	12:52 PM																		
256	Dilution	06/07/98	12:52 PM	06/07/98	01:22 PM																		
257	Agitate	06/07/98	01:22 PM	06/07/98	01:52 PM																		
258	CIP	06/07/98	01:52 PM	06/07/98	01:52 PM																		
259	SIP	06/07/98	01:52 PM	06/07/98	01:52 PM																		
260	Clean Up	06/07/98	01:52 PM	06/07/98	01:52 PM																		
261	Subtotal																						
262																							
263																							
264	12 A Cont. Cent./Solids																						
265	Set Up	06/07/98	12:52 PM	06/07/98	01:52 PM																		
266	Centrifugation	06/07/98	01:52 PM	06/07/98	02:22 PM																		
267	Wash	06/07/98	02:22 PM	06/07/98	02:28 PM																		
268	CIP	06/07/98	02:28 PM	06/07/98	02:28 PM																		
269	SIP	06/07/98	02:28 PM	06/07/98	02:28 PM																		
270	Clean Up	06/07/98	02:28 PM	06/07/98	02:28 PM																		
271	Sub Total																						
272																							
273																							
274	11 B Resolubilization																						
275	Set Up	06/07/98	02:28 PM	06/07/98	02:28 PM																		
276	Dilution	06/07/98	02:28 PM	06/07/98	02:58 PM																		
277	Agitate	06/07/98	02:58 PM	06/07/98	03:13 PM																		
278	CIP	06/07/98	03:13 PM	06/07/98	04:13 PM																		
279	SIP	06/07/98	04:13 PM	06/07/98	05:13 PM																		
280	Clean Up	06/07/98	05:13 PM	06/07/98	06:13 PM																		
281	Subtotal																						
282																							
283																							
284	12 B Cont. Cent./Solids																						
285	Set Up	06/07/98	02:13 PM	06/07/98	03:13 PM																		
286	Centrifugation	06/07/98	03:13 PM	06/07/98	03:43 PM																		
287	Wash	06/07/98	03:43 PM	06/07/98	03:49 PM																		
288																							

FIG. 45F

[illegible]

Fig. 456

QC Load Table - PE Module

	Operation	QA/QC Samples										Immunological	Act.											
		Start		Finish		Visual		Chemical		Biochemical														
		Date	Time	Date	Time	AV-1	AV-2	AC-1	AC-2	AC-3	AC-4			AC-5	AC-6	AB-1	AB-2	AB-3	AB-4	AB-5	AB-6	AB-7	AI-1	AI-2
338	Elite A	06/03/96	08:00 AM																					
339	Elite B	06/08/96	01:52 PM	06/08/96	03:12 PM																			
340	Regenerate	06/08/96	03:12 PM	06/08/96	03:12 PM																			
341	Store	06/08/96	03:12 PM	06/08/96	03:25 PM																			
342	CIP	06/08/96	03:25 PM	06/08/96	03:52 PM																			
343	SIP	06/08/96	04:52 PM	06/08/96	04:52 PM																			
344	Clean Up	06/08/96	05:52 PM	06/08/96	06:52 PM																			
345	Sub Total																							
346																								
347																								
348	17 A PIA MPLC																							
349	Equilibration	06/08/96	02:59 PM	06/08/96	03:38 PM																			
350	Load	06/08/96	03:12 PM	06/08/96	04:17 PM																			
351	Wash	06/08/96	04:17 PM	06/08/96	05:03 PM																			
352	Elite A	06/08/96	05:03 PM	06/08/96	05:49 PM																			
353	Elite B	06/08/96	05:49 PM	06/08/96	05:49 PM																			
354	Regenerate	06/08/96	05:49 PM	06/08/96	05:57 PM																			
355	Store	06/08/96	05:57 PM	06/08/96	06:13 PM																			
356	CIP	06/08/96	06:13 PM	06/08/96	07:13 PM																			
357	SIP	06/08/96	07:13 PM	06/08/96	08:13 PM																			
358	Clean Up	06/08/96	08:13 PM	06/08/96	08:13 PM																			
359	Sub Total																							
360																								
361																								
362	18 A Flow Dialysis																							
363	Set Up	06/08/96	03:29 PM	06/08/96	04:29 PM																			
364	Flush	06/08/96	04:29 PM	06/08/96	05:09 PM																			
365	Prime	06/08/96	05:09 PM	06/08/96	05:49 PM																			
366	Dialysis	06/08/96	05:49 PM	06/08/96	06:49 PM																			
367	Wash	06/08/96	06:49 PM	06/08/96	06:49 PM																			
368	Flush	06/08/96	06:49 PM	06/08/96	07:09 PM																			
369	Store	06/08/96	07:09 PM	06/08/96	07:49 PM																			
370	CIP	06/08/96	07:49 PM	06/08/96	08:49 PM																			
371	SIP	06/08/96	08:49 PM	06/08/96	09:49 PM																			
372	Clean Up	06/08/96	09:49 PM	06/08/96	10:49 PM																			
373	Sub Total																							
374																								
375																								
376	19 A PIA MPLC																							
377	Equilibration	06/08/96	05:59 PM	06/08/96	06:31 PM																			
378	Load	06/08/96	06:49 PM	06/08/96	07:03 PM																			
379	Wash	06/08/96	07:03 PM	06/08/96	07:41 PM																			
380	Elite A	06/08/96	07:41 PM	06/08/96	08:20 PM																			
381	Elite B	06/08/96	08:20 PM	06/08/96	08:20 PM																			
382	Regenerate	06/08/96	08:20 PM	06/08/96	08:26 PM																			
383	Store	06/08/96	08:26 PM	06/08/96	08:39 PM																			
384	CIP	06/08/96	08:39 PM	06/08/96	09:39 PM																			
385	SIP	06/08/96	09:39 PM	06/08/96	10:39 PM																			
386	Sub Total																							

FIG. 4SH

QC Load Table - PE Module

QC Load Table - PE Module

	Operation	Start				Finish		QA/QC Samples												Immunochemical				Act.			
		Date		Time		Date		Time		Visual			Chemical			Biochemical						AB-1			AB-2		
										AV-1	AV-2	AC-1	AC-2	AC-3	AC-4	AC-5	AC-6	AB-1	AB-2	AB-3	AB-4	AB-5	AB-6	AB-7	AI-1	AI-2	AA-1
387	Clean Up	06/03/96	08:00 AM	06/08/96	10:39 PM	06/08/96	11:39 PM																				
388	Sub Total																										
389																											
390	20 A Flow Dialysis																										
391	Set Up	06/08/96	07:00 PM	06/08/96	07:00 PM	06/08/96	07:00 PM																				
392	Flush	06/08/96	07:00 PM	06/08/96	07:40 PM	06/08/96	08:20 PM																				
393	Prime	06/08/96	07:40 PM	06/08/96	08:20 PM	06/08/96	10:20 PM																				
394	Dialysis	06/08/96	08:20 PM	06/08/96	10:20 PM	06/08/96	10:20 PM																				
395	Wash	06/08/96	10:20 PM	06/08/96	10:20 PM	06/08/96	10:20 PM																				
396	Flush	06/08/96	10:20 PM	06/08/96	10:40 PM	06/08/96	11:20 PM																				
397	Store	06/08/96	10:40 PM	06/08/96	11:20 PM	06/08/96	11:20 PM																				
398	CIP	06/08/96	11:20 PM	06/08/96	11:20 PM	06/08/96	11:20 PM																				
399	SIP	06/08/96	11:20 PM	06/08/96	11:20 PM	06/08/96	11:20 PM																				
400	Clean Up	06/08/96	11:20 PM	06/08/96	12:20 AM	06/08/96	12:20 AM																				
401	Sub Total																										
402																											
403	21 A PIA MPLC																										
404	Equilibration	06/08/96	09:28 PM	06/08/96	09:57 PM	06/08/96	10:28 PM																				
405	Load	06/08/96	10:20 PM	06/08/96	10:26 PM	06/08/96	11:01 PM																				
406	Wash	06/08/96	10:26 PM	06/08/96	11:01 PM	06/08/96	11:36 PM																				
407	Elute A	06/08/96	11:01 PM	06/08/96	11:36 PM	06/08/96	11:36 PM																				
408	Elute B	06/08/96	11:36 PM	06/08/96	11:36 PM	06/08/96	11:42 PM																				
409	Regenerate	06/08/96	11:36 PM	06/08/96	11:42 PM	06/08/96	11:54 PM																				
410	Store	06/08/96	11:42 PM	06/08/96	11:54 PM	06/08/96	11:54 PM																				
411	CIP	06/08/96	11:54 PM	06/08/96	11:54 PM	06/08/96	11:54 PM																				
412	SIP	06/08/96	11:54 PM	06/08/96	11:54 PM	06/08/96	11:54 PM																				
413	Clean Up	06/08/96	11:54 PM	06/08/96	12:54 AM	06/08/96	12:54 AM																				
414	Sub Total																										
415																											
416	22 A Sterile Filtration																										
417	Set Up	06/09/96	08:06 AM	06/09/96	08:36 AM	06/09/96	12:08 AM																				
418	Filtration	06/08/96	11:36 PM	06/09/96	12:06 AM	06/09/96	12:36 AM																				
419	Storage	06/09/96	12:06 AM	06/09/96	12:36 AM	06/09/96	12:36 AM																				
420	CIP	06/09/96	12:36 AM	06/09/96	12:36 AM	06/09/96	12:36 AM																				
421	SIP	06/09/96	12:36 AM	06/09/96	12:36 AM	06/09/96	12:36 AM																				
422	Clean Up	06/09/96	12:36 AM	06/09/96	01:36 AM	06/09/96	01:36 AM																				
423	Sub Total																										
424																											
425																											
426																											

FIG. 45I




FIG. 46




FIG. 47




FIG. 48




FIG. 49












FIG. 50









```
graph TD; 4706[/4706  
prep equipment list  
with functional  
specifications/] --> 5904[generate prep equipment  
calibration tables]; 5902[/5902  
prep equipment  
calibration  
data/] --> 5904; 5904 --> 5906[/5906  
prep equipment  
calibration  
tables/];
```

FIG. 59

FIG. 59




FIG. 60




FIG-61




FIG-62

Equipment Maintenance Table - Microbial Fermentation

6408

6406

6404

6402

Equipment Items	Fillers						Gaskets						Bearings					
	Materials			Labor			Materials			Labor			Materials					
	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.			
Inoculum Freezer																		
-80 C Stock Freezer																		
Shaking Water Bath																		
Flow Incubator-Shaker																		
Microscope																		
Seed Fermentation																		
Seed Bioreactor																		
Production Bioreactor	75868	1	100	55	.55	.5	.0875	48994	1	500	55	.11	1	.035				
Harvest Heat Exchanger								62589	1	350	85	.2429	1	.05				
Harvest Vessel																		
Agitator																		
Filter Holder																		
Manifolding																		
Instrumentation																		
MF Flush Vessel																		
MF Prime Vessel																		
MF Filtrate Vessel																		
Agitator																		
MF Wash Vessel																		
MF Regeneration Vessel																		
MF Storage Vessel																		

FIG. 64A

Equipment Maintenance Table - Microbial Fermentation

6408

6410

6412

Equipment Items	Seals						Belts					
	Labor			Materials			Labor			Materials		
	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle
Production Bioreactor												
Harvest Heat Exchanger												
Harvest Vessel												
Agitator												
Pump												
Filter Holder												
Manifolding												
Instrumentation												
MF Flush Vessel												
MF Prime Vessel												
MF Filtrate Vessel												
Agitator												
MF Wash Vessel												
MF Regeneration Vessel												
MF Storage Vessel												
Seed Bioreactor												
Microscope												
Floor Incubator-Shaker												
Shaking Water Bath												
-80 C Stock Freezer												
Production Bioreactor												
Harvest Heat Exchanger												
Harvest Vessel												
Agitator												
Pump												
Filter Holder												
Manifolding												
Instrumentation												
MF Flush Vessel												
MF Prime Vessel												
MF Filtrate Vessel												
Agitator												
MF Wash Vessel												
MF Regeneration Vessel												
MF Storage Vessel												

FIG. 64B

Equipment Maintenance Table - Microbial Fermentation

6418

6416

Equipment Items	Shafts										Lubricant					
	Labor					Materials					Labor			Materials		
	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	
Inoculum Prep																
-80 C Stock Freezer																
Shaking Water Bath																
22 flask growing																
Floor Incubator-Shaker																
Microscope																
23 Seed Fermentation																
Seed Bioreactor													78154	.5		
24 Fermentation	500	25	.05	1	.035											
Production Bioreactor																
25 Whole Cell Harvest																
Harvest Heat Exchanger																
Harvest Vessel																
Agitator																
26 Cell Concentration																
Pump																
Filter Holder																
Manifolding																
Instrumentation																
MF Flush Vessel																
MF Prime Vessel																
MF Filtrate Vessel																
Agitator																
MF Wash Vessel																
MF Regeneration Vessel																
MF Storage Vessel																

Fig. 64C

Equipment Maintenance Table - Microbial Fermentation

6418

6420

Equipment Items	Thermal Media									
	Labor			Materials			Labor			
	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours
Inoculum Prep										
-80 C Stock Freezer										
Shaking Water Bath										
Flow Incubator-Shaker										
Microscope										
Seed Bioreactor	1.5	.03	.5	.175						
Production Bioreactor					56258	5	500	.85	425	1
Whole Cell Harvest										
Harvest Heat Exchanger										
Harvest Vessel										
Agitator										
Cell Concentration										
Pump										
Filter Holder										
Manifolding										
Instrumentation										
MF Flush Vessel										
MF Prime Vessel										
MF Filtrate Vessel										
Agitator										
MF Wash Vessel										
MF Regeneration Vessel										
MF Storage Vessel										

Fig. 64D

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Filters						Gaskets						Bearings			
	Materials			Labor			Materials			Labor			Materials			
	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	Item No.	\$/Cycle	Hours	Item No.
MF Wash Vessel																
Pump																
Filter Holder																
Manifolding																
Instrumentation																
MF Flush Vessel																
MF Prime Vessel																
MF Filtrate Vessel																
MF Wash Vessel																
MF Regeneration Vessel																
MF Storage Vessel																
Resuspension Vessel																
Stir Plate																
Cell Disruptor																
Lysate Vessel																
Resuspension Vessel																
Stir Plate																
MF Wash Vessel																
Pump																
Filter Holder																

FIG-64E

Equipment Maintenance Table Microbial Fermentation

Equipment Items	Seals						Belts					
	Labor			Materials			Labor			Materials		
	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle
Cell Concentration												
MF Wash Vessel												
Pump												
Filter Holder												
Manifolding												
Instrumentation												
MF Flush Vessel												
MF Prime Vessel												
MF Filtrate Vessel												
MF Wash Vessel												
MF Regeneration Vessel												
MF Storage Vessel												
Cell Resuspension												
Resuspension Vessel												
Stir Plate												
Cell Disruptor												
Lysate Vessel												
Cell Resuspension												
Resuspension Vessel												
Stir Plate												
Cell Concentration												
MF Wash Vessel												
Pump												
Filter Holder												

Fig. 64F

Equipment Maintenance Table for Microbial Fermentation

Equipment Items	Shafts						Lubricant					
	Labor			Materials			Labor			Materials		
	Unit Cost	\$/Cycle	Hours	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	Item No.	Qty	Cycle Life
MF Wash Vessel												
Pump												
Filter Holder												
Manifolding												
Instrumentation												
MF Flush Vessel												
MF Prime Vessel												
MF Filtrate Vessel												
MF Wash Vessel												
MF Regeneration Vessel												
MF Storage Vessel												
Resuspension Vessel												
Stir Plate												
Cell Disruptor												
Lysate Vessel												
Resuspension Vessel												
Stir Plate												
MF Wash Vessel												
Pump												
Filter Holder												

FIG. 64G

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Filters						Gaskets						Bearings	
	Materials			Labor			Materials			Labor			Materials	
	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	Item No.	
Manifolding														
Instrumentation														
MF Flush Vessel														
MF Prime Vessel														
MF Filtrate Vessel														
MF Dilute Vessel														
MF Wash Vessel														
MF Regeneration Vessel														
MF Storage Vessel														
UF Renaturation Vessel														
Renaturation Vessel														
Stir Plate														
Buffer Exchange Vessel														
Pump														
Filter Holder														
Manifolding														
Instrumentation														
UF Flush Vessel														
UF Prime Vessel														
UF Filtrate Vessel														
UF Wash Vessel														
UF Diluent Vessel														
UF Regeneration Vessel														
UF Storage Vessel														

FIG-64I

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Seals						Bells			
	Labor			Materials			Labor		Materials	
	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Qty	Cycle Life	Unit Cost	\$/Cycle
Manifolding										
Instrumentation										
MF Flush Vessel										
MF Prime Vessel										
MF Filtrate Vessel										
MF Dilute Vessel										
MF Wash Vessel										
MF Regeneration Vessel										
MF Storage Vessel										
Renaturant Vessel										
Stir Plate										
Pump										
Filter Holder										
Manifolding										
Instrumentation										
UF Flush Vessel										
UF Prime Vessel										
UF Filtrate Vessel										
UF Wash Vessel										
UF Diluent Vessel										
UF Regeneration Vessel										
UF Storage Vessel										

FIG-64J

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Shafts						Lubricant			
	Labor			Materials			Materials			
	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours
	Cycle Life									
Manifolding										
Instrumentation										
MF Flush Vessel										
MF Prime Vessel										
MF Filtrate Vessel										
MF Dilute Vessel										
MF Wash Vessel										
MF Regeneration Vessel										
MF Storage Vessel										
Renaturation										
Renaturation Vessel										
Stir Plate										
Buffer Exchange										
Pump										
Filter Holder										
Manifolding										
Instrumentation										
UF Flush Vessel										
UF Prime Vessel										
UF Filtrate Vessel										
UF Wash Vessel										
UF Diluent Vessel										
UF Regeneration Vessel										
UF Storage Vessel										

Fig. 64 K

Equipment Maintenance Table- Microbial Fermentation

Equipment Items	Thermal Media							
	Labor				Materials			
	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost
Manifolding								
Instrumentation								
MF Flush Vessel								
MF Prime Vessel								
MF Filtrate Vessel								
MF Dilute Vessel								
MF Wash Vessel								
MF Regeneration Vessel								
MF Storage Vessel								
Renaturation Vessel								
Renaturation Vessel								
Stir Plate								
Pump								
Filter Holder								
Manifolding								
Instrumentation								
UF Flush Vessel								
UF Prime Vessel								
UF Filtrate Vessel								
UF Wash Vessel								
UF Diluent Vessel								
UF Regeneration Vessel								
UF Storage Vessel								

FIG. 64L

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Filters						Gaskets						Bearings		
	Materials			Labor			Materials			Labor			Materials		
	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	Item No.	\$/Cycle	Item No.
UF Waste Vessel															
Chromatography Column															
Pump															
Inst. & Control System															
Manifolding															
Equilibration Vessel															
Wash Vessel															
Eluent Vessel															
Regenerate Vessel															
Storage Vessel															
Waste Vessel (1)															
Product Vessel															
Waste Vessel (2)															
Chromatography Column															
Pump															
Inst. & Control System															
Manifolding															
Equilibration Vessel															
Wash Vessel															
Eluent Vessel															
Regenerate Vessel															

FIG. 64h

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Seals								Belts							
	Labor				Materials				Labor				Materials			
	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	
UF Waste Vessel																
Chromatography Column																
Pump																
Inst. & Control System																
Manifolding																
Equilibration Vessel																
Wash Vessel																
Eluent Vessel																
Regenerate Vessel																
Storage Vessel																
Waste Vessel (1)																
Product Vessel																
Waste Vessel (2)																
Chromatography Column																
Pump																
Inst. & Control System																
Manifolding																
Equilibration Vessel																
Wash Vessel																
Eluent Vessel																
Regenerate Vessel																

FIG. 64N

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Shafts						Lubricant								
	Labor			Materials			Labor			Materials					
	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life
UF Waste Vessel															
Chromatography Column															
Pump															
Inst. & Control System															
Manifolding															
Equilibration Vessel															
Wash Vessel															
Eluent Vessel															
Regenerate Vessel															
Storage Vessel															
Waste Vessel (1)															
Product Vessel															
Waste Vessel (2)															
Chromatography Column															
Pump															
Inst. & Control System															
Manifolding															
Equilibration Vessel															
Wash Vessel															
Eluent Vessel															
Regenerate Vessel															

FIG. 640

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Thermal Media										
	Labor			Materials							
	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle
UF Waste Vessel											
Chromatography Column											
Pump											
Inst. & Control System											
Manifolding											
Equilibration Vessel											
Wash Vessel											
Eluent Vessel											
Regenerate Vessel											
Storage Vessel											
Waste Vessel (1)											
Product Vessel											
Waste Vessel (2)											
Chromatography Column											
Pump											
Inst. & Control System											
Manifolding											
Equilibration Vessel											
Wash Vessel											
Eluent Vessel											
Regenerate Vessel											

FIG. 64P

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Shafts								Lubricant		
	Labor				Materials				Labor		
	Unit Cost	\$/Cycle	Hours	Cycle Life	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	Cycle Life
Storage Vessel											
Waste Vessel (1)											
Product Vessel											
Waste Vessel (2)											
Buffer Exchange											
Pump											
Filter Holder											
Manifolding											
Instrumentation											
UF Flush Vessel											
UF Prime Vessel											
UF Filtrate Vessel											
UF Wash Vessel											
UF Diluent Vessel											
UF Regeneration Vessel											
UF Storage Vessel											
UF Waste Vessel											
Chromatography Column											
Pump											
Inst. & Control System											
Manifolding											
Equilibration Vessel											

Fig. 64

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Thermal Media									
	Labor					Materials				
	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours
Storage Vessel										
Waste Vessel (1)										
Product Vessel										
Waste Vessel (2)										
Buffer Exchange										
Pump										
Filter Holder										
Manifolding										
Instrumentation										
UF Flush Vessel										
UF Prime Vessel										
UF Filtrate Vessel										
UF Wash Vessel										
UF Diluent Vessel										
UF Regeneration Vessel										
UF Storage Vessel										
UF Waste Vessel										
Chromatography Column										
Pump										
Inst. & Control System										
Manifolding										
Equilibration Vessel										

FIG. 64T

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Filters							Gaskets							Bearings		
	Materials				Labor			Materials				Labor			Materials		Item No.
	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle			
Wash Vessel																	
Eluent Vessel																	
Regenerate Vessel																	
Storage Vessel																	
Waste Vessel (1)																	
Product Vessel																	
Waste Vessel (2)																	
200 Buffer Exchange																	
Pump																	
Filter Holder																	
Manifolding																	
Instrumentation																	
UF Flush Vessel																	
UF Prime Vessel																	
UF Filtrate Vessel																	
UF Wash Vessel																	
UF Diluent Vessel																	
UF Regeneration Vessel																	
UF Storage Vessel																	
UF Waste Vessel																	
200 Chromatography																	
Chromatography Column																	
Pump																	

Fig. 64U

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Seals						Belts								
	Labor			Materials			Labor			Materials					
	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty
Wash Vessel															
Eluent Vessel															
Regenerate Vessel															
Storage Vessel															
Waste Vessel (1)															
Product Vessel															
Waste Vessel (2)															
200 Buffer Exchange															
Pump															
Filter Holder															
Manifolding															
Instrumentation															
UF Flush Vessel															
UF Prime Vessel															
UF Filtrate Vessel															
UF Wash Vessel															
UF Diluent Vessel															
UF Regeneration Vessel															
UF Storage Vessel															
UF Waste Vessel															
21 Chromatography															
Chromatography Column															
Pump															

FIG. 64V

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Shafts						Lubricant								
	Labor			Materials			Labor			Materials					
	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life
Wash Vessel															
Eluent Vessel															
Regenerate Vessel															
Storage Vessel															
Waste Vessel (1)															
Product Vessel															
Waste Vessel (2)															
200 Buffer Exchange															
Pump															
Filter Holder															
Manifolding															
Instrumentation															
UF Flush Vessel															
UF Prime Vessel															
UF Filtrate Vessel															
UF Wash Vessel															
UF Diluent Vessel															
UF Regeneration Vessel															
UF Storage Vessel															
UF Waste Vessel															
200 Chromatography															
Chromatography Column															
Pump															

FIG. 64W

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Thermal Media									
	Labor			Materials						
	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Labor
Wash Vessel										
Eluent Vessel										
Regenerate Vessel										
Storage Vessel										
Waste Vessel (1)										
Product Vessel										
Waste Vessel (2)										
20 Buffer Exchange										
Pump										
Filter Holder										
Manifolding										
Instrumentation										
UF Flush Vessel										
UF Prime Vessel										
UF Filtrate Vessel										
UF Wash Vessel										
UF Diluent Vessel										
UF Regeneration Vessel										
UF Storage Vessel										
UF Waste Vessel										
21 Chromatography										
Chromatography Column										
Pump										

Fig. 64 X

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Filters						Gaskeits						Bearings		
	Materials			Labor			Materials			Labor			Materials		
	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.
Inst. & Control System															
Manifolding															
Equilibration Vessel															
Wash Vessel															
Eluent Vessel															
Regenerate Vessel															
Storage Vessel															
Waste Vessel (1)															
Product Vessel															
Waste Vessel (2)															
Station Filtration															
MF Wash Vessel															
Pump															
Filter Holder															
Manifolding															
Instrumentation															
MF Flush Vessel															
MF Prime Vessel															
MF Filtrate Vessel															
MF Wash Vessel															

00403000 : 004000

FIG. 64Y

[illegible]

6/27/2017

Equipment Maintenance Table - Microbial Fermentation

[illegible]

1. The first step is to identify the problem or question that needs to be addressed. This involves understanding the context and the specific requirements of the task.

FI 64 AA

Equipment Maintenance Table - Microbial Fermentation

Equipment Items	Thermal Media									
	Labor				Materials					
	Unit Cost	\$/Cycle	Hours	\$/Cycle	Item No.	Qty	Cycle Life	Unit Cost	\$/Cycle	Hours
Inst. & Control System										
Manifolding										
Equilibration Vessel										
Wash Vessel										
Eluent Vessel										
Regenerate Vessel										
Storage Vessel										
Waste Vessel (1)										
Product Vessel										
Waste Vessel (2)										
MF Wash Vessel										
Pump										
Filter Holder										
Manifolding										
Instrumentation										
MF Flush Vessel										
MF Prime Vessel										
MF Filtrate Vessel										
MF Wash Vessel										

FIG. 64A-B

Master Process Parameters Table - Biopharmaceutical

Unit Operation Type	Group 1			Group 2			Group 3		
	Parameter	Soln.		Parameter	Soln.		Parameter	Soln.	
T1 Inoculum Prep	Number of Flasks Media Volume/Flask	2 0.25 Liters		Temperature Agitation Duration	37 C 200 RPM 18 Hours		Final OD		12
T2 Flask Growth	Scale Up Ratio Media Volume/Flask	10 Fold 1.25 L		Temperature Agitation Duration	37 C 200 RPM 18 Hours		Final OD		12
T3 Fermentation Production	Scale Up Ratio Fermentor Working Volume Antifoam A Antifoam B Base Add	S-101 500 Liters 1 M/L 1 M/L 5 M/L 5 M/L		Growth Temperature Agitation Spurge Rate Back Pressure Total Duration	37 Hours 1 HP/100L 1.5 VVM 5 PSIG 21 Hrs		Final OD Dry Cell Mass Product Concentration CIP		12 9.86 Gm TDCM/L 0.3 Gms Product/L Y
T4 Initial seeding	Number of Ampules Volume Per Ampule Saring Cell Density Ampule Split Ratio Culture Vessel Type Feed Volume	2 300,000 Cells/Ml 1 Vessels/Ampule Roll. Bot. 100 Ml		Serum Content Feed Rate Days to Confluence	2.0% Fetal Bovine Serum 1 Feed per vessel per 2 Days 2 Days		Amplification Factor		100%
T5 Culture Vessel Split	Vessel Split Ratio New Vessel Type Feed Volume Serum Content	2 100 Ml 2.0% Fetal Bovine Serum		Feed Rate Days to Confluence	1 Feed per vessel per 2 Days 2 Days		Amplification Factor		100%
T6 Spinner Flask Seeding	Flask Feed Volume Vessel/Flask Ratio uCarrier Density Number of PBS Washes Number of Media Washes No. of Media/Serum Washes	4 Liters 0.1 L Cells/L Flask 5 Gm/Liter 2 1 2 FBS		Serum Content Feed Rate Days to Confluence	2.0% Fetal Bovine Serum 1 Feed per vessel per 2 Days 2 Days		Amplification Factor		100%
T7 Biocytinthesis Bioreactor Preparation (Stirred Tank Reactor)	Reactor Feed Volume Spinner/Reactor Ratio uCarrier Density Number of PBS Washes Number of Media Washes No. of Media/Serum Washes	500 Liters 8.3 5 Gm/Liter 2 1 2		Serum Content Feed Rate Days to Confluence Serum Free Media Washes	2.0% Fetal Bovine Serum 1 Feed per vessel per 2 Days 10 Days 2		Product Concentration Total Protein Concn.		2500% Mg Prod/L 0.125 Mg TPMI
T8 Biocytinthesis Bioreactor Preparation (Hollow Fiber Reactor)	Reactor Feed Volume Number of PBS Washes Number of Media Washes No. of Media/Serum Washes Serum Content	100 Liters 2 2 2 2.0% Fetal Bovine Serum		Number of Reactors Feed Rate Days to Confluence	1 Feed per vessel per 1 Days 10 Days		Harvest Volume Product Concentration Total Protein Concn.		500% Liters 25 Mg Prod/L 0.125 Mg TPMI
T9 Biocytinthesis Bioreactor Preparation (Fluidized Bed Reactor)	Reactor Feed Volume uCarrier Density Number of PBS Washes Number of Media Washes No. of Media/Serum Washes Serum Content	Liters Gms/L		Number of Reactors Feed Rate Days to Confluence	1 Feed per vessel per 1 Days 10 Days		Product Concentration Total Protein Concn.		2500% Mg Prod/L 0.125 Mg TPMI
T10 Initial seeding	Number of Ampules Volume Per Ampule Saring Cell Density Ampule Split Ratio	2 300,000 Cells/Ml 1 Vessels/Ampule		Serum Content Feed Rate Days to Confluence	2.0% Fetal Bovine Serum 1 Feed per vessel per 2 Days 2 Days		Amplification Factor		100%

Master Process Parameters Table - Biopharmaceutical

Unit Operation Type	Group 1		Group 2		Group 3	
	Parameter	Soln.	Parameter	Soln.	Parameter	Soln.
	Culture Vessel Type		PBS Washes			
	Feed Volume	Rot. Bot. 100 MI	Trypan Wash			
	Vessel Spill Ratio	2	Feed Rate			
	New Vessel Type	RB	Days to Confluence			
	Feed Volume	100 MI	PBS Washes			
	Serum Content	2.0% Fetal Bovine Serum	Trypan Wash			
T11 Culture Vessel Spill					Amplification Factor	100%
	Flask Feed Volume	4 Liters	Serum Content			
	Vessel/Flask Ratio	0.1 L Cells/L Flask	Feed Rate			
	Carrier Density	5 Gm/Liter	Days to Confluence			
	Number of PBS Washes	2				
	Number of Media Washes	1				
	No. of Media/Serum Washes	2				
T12 Spinner Flask Spill					Amplification Factor	100%
	Reactor Feed Volume	500 Liters	Serum Content			
	Spinner/Reactor Ratio	8.3	Feed Rate			
	Carrier Density	5 Gm/Liter	Days to Confluence			
	Number of PBS Washes	2				
	Number of Media Washes	1				
	No. of Media/Serum Washes	2				
T13 Biosynthesis Bioreactor Preparation (Stirred Tank Reactor)					Product Concentration Total Protein Conc.	2500% Mg Prod/L 0.125 Mg TP/MI
	Reactor Feed Volume	500 Liters	Serum Content			
	Spinner/Reactor Ratio	8.3	Feed Rate			
	Carrier Density	5 Gm/Liter	Days to Confluence			
	Number of PBS Washes	2				
	Number of Media Washes	1				
	No. of Media/Serum Washes	2				
T14 Biosynthesis Bioreactor Preparation (Fluidized Bed Reactor)					Product Concentration Total Protein Conc.	2500% Mg Prod/L 0.125 Mg TP/MI
	Reactor Feed Volume	500 Liters	Serum Content			
	Spinner/Reactor Ratio	8.3	Feed Rate			
	Carrier Density	5 Gm/Liter	Days to Confluence			
	Number of PBS Washes	2				
	Number of Media Washes	1				
	No. of Media/Serum Washes	2				
T15 Initial Coupling					Amplification Factor	100%
	Flask Feed Volume	4 Liters	Serum Content			
	Vessel/Flask Ratio	0.1 L Cells/L Flask	Feed Rate			
	Carrier Density	5 Gm/Liter	Days to Confluence			
	Number of PBS Washes	2				
	Number of Media Washes	1				
	No. of Media/Serum Washes	2				
T16 Additional Coupling					Product Concentration Total Protein Conc.	2500% Mg Prod/L 0.125 Mg TP/MI
	Reactor Feed Volume	500 Liters	Serum Content			
	Spinner/Reactor Ratio	8.3	Feed Rate			
	Carrier Density	5 Gm/Liter	Days to Confluence			
	Number of PBS Washes	2				
	Number of Media Washes	1				
	No. of Media/Serum Washes	2				
T17 Peptide Cleavage					Harvest Volume Product Concentration Total Protein Conc.	500% Liters 25 Mg Prod/L 0.125 Mg TP/MI
	Reactor Feed Volume	100 Liters	Serum Content			
	Spinner/Reactor Ratio	2	Feed Rate			
	Carrier Density	2	Days to Confluence			
	Number of PBS Washes	2				
	Number of Media Washes	1				
	No. of Media/Serum Washes	2				
T18 Tissue Thawing					Temperature Regulation CIP SIP	Y Y Y
	Crude Product Yield	25 Gm Crude Prod/Kg Tissue	Contaminant Protein Conc.			
	Environmental Temperature	25 C				
	Thaw Duration	16 Hours				
T19 Homogenization					Temperature Regulation CIP SIP	Y Y Y
	Crude Product Yield	25 Gm Crude Prod/Kg Tissue	Contaminant Protein Conc.			
	Liquid/Solid Ratio	10 L Solution/Kg Tissue				
	Homogenization Temp.	4 C				
	Homogenizer Type	RS				
	Energy Input	200 HP/100L/Hr				
	Duration	4 Hours				
T20 Liquid Thawing					Amplification Factor	100%

Master Process Parameters Table - Biopharmaceutical

Unit Operation Type	Group 1		Group 2		Group 3	
	Parameter	Soln.	Parameter	Soln.	Parameter	Soln.
T21	Product Ppt by Solids		Reagent Concentration	1 M	Kgms of Reagent/Liters Product Temperature Addition Time Additional Mix Time	0.25 Kg/L 4 C 0.5 Hours 2 Hours
					Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	95% 95% Y Y
T22	Product Ppt by Liquids		Reagent Concentration	1 M	Liters Reagent/Liters Product Temperature Addition Time Additional Mix Time	0.25 L/L 4 C 0.5 Hours 2 Hours
					Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	95% 95% Y Y
T23	Contaminant Ppt by Solids		Reagent Concentration	1 M	Kgms of Reagent/Liters Product Temperature Addition Time Additional Mix Time	0.25 Kg/L 4 C 0.5 Hours 2 Hours
					Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	95% 95% Y Y
T24	Contaminant Ppt by Liquids		Reagent Concentration	1 M	Liters Reagent/Liters Product Temperature Addition Time Additional Mix Time	0.25 L/L 4 C 0.5 Hours 2 Hours
					Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	95% 95% Y Y
T25	Solids Harvest Tangential Flow MF		Porosity Average Flux Rate Concentration Factor Wash Regenerate Store	0.2 Micron 11 L/SF/HR at 40 Psg at 4 C 400 Liters/SF 1 HR	Flush Prime Wash Regenerate Store	2 L/SF 2 L/SF 10 Fold 0.5 L/SF 1 L/SF 2 L/SF
					Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	95% 95% Y Y
T26	Continuous Centrifugation Solids Harvest		System Void Volume	5 Liters	RCF Time Volume Reduction Wash Volume	10,000 X G 60 Minutes 30 X Vol. Reduction 0.2 X System Void Volume
					Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	95% 95% Y Y
T27	Continuous Centrifugation Supernatant Harvest		System Void Volume	6 Liters	RCF Time Volume Reduction Wash Volume	10,000 X G 30 Minutes 0.002 Vol. Reduction 1.5 X System Void Volume
					Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	95% 0.3 Y Y
T28	Dilution		System Void Volume	6 Liters	RCF Time Volume Reduction Wash Volume	10,000 X G 30 Minutes 18 X Vol. Reduction 1.5 X System Void Volume
					Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	95% 0.95 Y Y
T29	Batch Centrifugation Solids Harvest		System Void Volume	6 Liters	RCF Time Volume Reduction Wash Volume	10,000 X G 30 Minutes 18 X Vol. Reduction 1.5 X System Void Volume
					Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP	95% 0.95 Y Y

Master Process Parameters Table - Biopharmaceutical

Unit Operation Type	Group 1			Group 2			Group 3		
	Parameter	Soln.		Parameter	Soln.		Parameter	Soln.	
T30 Batch Centrifugation Supernatant Harvest	System Void Volume		6 Liters	RCF Time Volume Reduction Wash Volume		10000 X G 30 Minutes 18 X Vol. Reduction 1.5 X System Void Volume	SIP Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP		Y 85% 95% Y Y Y
T31 Cell Disruption High Press. Homogen.	Product Temperature Utility Temperature Void Volume		8 C 2 C 5 Liters	Number of Passes Pressure Flow Rate Temperature Increase		6 Times 12,000 PSI 5 LPM 1.8 Degrees C/1,000 PSI	Rinse Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP		500% Void Volumes 95% 95% Y Y Y
T32 Cell Disruption Bead Mill	Number of Passes Bead Size Void Volume Flow Rate		2 0.5 LPM				Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP		95% Y Y Y Y
T33 Cell Disruption Chemical Lysis	Reagent Temperature Exposure Time		0.5 M NaOH 4 C 2 Hours	Urea Reagent/Gm Product Titration		0.4 L/Gm 0 M/Liter	Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP		95% Y Y Y Y
T34 Microfiltration Tangential Flow	Porosity Average Flux Rate Total Throughput Filtration Time		0.2 Micron 50 LSF/HR at 40 Palg at 4 C 400 L/hr/SF 2 HR	Flush Prime Wash Solids Regenerate Store		2.00 LSF 2.00 LSF 0.50 LSF 0.30% Of Product Solution 1.00 LSF 2.00 LSF	Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP		95% 95% Y Y Y Y
T35 Microfiltration Dried End	Porosity Average Flux Rate Total Throughput Filtration Time		0.2 Micron 50 LSF/HR at 40 Palg at 4 C 400 L/hr/SF 0.5 HR	Flush Prime Wash Solids Regenerate Store		0 LSF 0 LSF 0.5 LSF 0.003 Of Product Solution 1 LSF 2 LSF	Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP		95% 0.95 N N N N
T36 Ultrafiltration Concentration/Clusion	Porosity Average Flux Rate Concentration Time		60 K MWL 3 LSF/HR at 40 Palg at 4 C 2 HR	Flush Prime Wash Dilute Concentrate Solids Regenerate		2.00 LSF 2.00 LSF 0.50 LSF 10.0 Fold 0.30% Of Product Solution 1.00 LSF	Store Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP		2.00 LSF 95% 95% Y Y Y
T37 Ultrafiltration Flow Dialysis	Porosity Average Flux Rate Dialysis Time		60 K MWL 3 LSF/HR at 40 Palg at 4 C 2 HR	Flush Prime Dialysis Buffer Wash Solids Regenerate		2 LSF 2.00 LSF 5.0 X Feed Stream Volume 0.50 LSF 0.30% Of Product Solution 1.00 LSF	Store Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP		200% LSF 95% 95% Y Y Y
T38 Prod. Ads. Chromatography HPLC	Column Capacity Column Oversize Factor Column Aspect Ratio Max. Linear Velocity		10 MG Prod./ML Of Packing 1.5 Fold 0.37 HD 100 Cm/Hr at 45 Palg and 4 C	Column Equilibration Column Wash Column Elute A Column Elute B Column Regenerate Column Store		5 Column Volumes 3 Column Volumes 3 Column Volumes 0 Column Volumes 1 Column Volumes 2 Column Volumes	Prod. Elution Volume Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP		80% 95% 95% N Y Y
T39 Prod. Ads. Chromatography MPLC	Column Capacity Column Oversize Factor		10 MG Prod./ML Of Packing 1.5 Fold	Column Equilibration Column Wash		5 Column Volumes 3 Column Volumes	Prod. Elution Volume Step Recovery of Product		80% 95%

A4

Master Process Parameters Table - Biopharmaceutical

Unit Operation Type	Group 1		Group 2		Group 3	
	Parameter	Soln.	Parameter	Soln.	Parameter	Soln.
	Column Aspect Ratio Max. Linear Velocity	0.37 HD 100 Cn/Hr at 45 Pa/g and 4 C	Column Elute A Column Elute B Column Regenerate Column Store		Step Recovery of T.P. Temperature Regulation CIP SIP	N Y Y Y
T140 Prod. Ads. Chromatography LPLC	Column Capacity Column Oversize Factor Column Aspect Ratio Max. Linear Velocity	10 MG Prod./MI Of Packing 1.5 Fold 0.37 HD 100 Cn/Hr at 45 Pa/g and 4 C	Column Equilibration Column Wash Column Elute A Column Elute B Column Regenerate Column Store		Prod. Elution Volume Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	42% 95% 95% N Y Y Y
T141 Cont. Ads. Chromatography HPLC	Column Capacity Column Oversize Factor Column Aspect Ratio Max. Linear Velocity	30 MG Cont./MI Of Packing 1.5 Fold 0.37 HD 100 Cn/Hr at 45 Pa/g and 4 C	Column Equilibration Column Wash Column Elute A Column Elute B Column Regenerate Column Store		Prod. Elution Volume Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	42% 95% 95% N Y Y Y
T142 Cont. Ads. Chromatography MPLC	Column Capacity Column Oversize Factor Column Aspect Ratio Max. Linear Velocity	10 MG Cont./MI Of Packing 1.5 Fold 0.37 HD 100 Cn/Hr at 45 Pa/g and 400% C	Column Equilibration Column Wash Column Elute A Column Elute B Column Regenerate Column Store		Prod. Elution Volume Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	42% 95% 95% N Y Y Y
T143 Cont. Ads. Chromatography LPLC	Column Capacity Column Oversize Factor Column Aspect Ratio Max. Linear Velocity	10 MG Cont./MI Of Packing 1.5 Fold 0.37 HD 100 Cn/Hr at 45 Pa/g and 4 C	Column Equilibration Column Wash Column Elute A Column Elute B Column Regenerate Column Store		Prod. Elution Volume Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	42% Columns Volumes 95% 95% N Y Y Y
T144 Size Excl. Chromatography HPLC	Load Capacity Length Max. Linear Velocity Void Volume	5% of Total Column Volume 100 Cn 100 Cn/Hr at 45 Pa/g and 4 C 25% Column Volume	Column Equilibration Column Wash Column Regenerate Column Store		Prod. Elution Volume Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	42% Columns Volumes 95% 95% N Y Y Y
T145 Size Excl. Chromatography MPLC	Load Capacity Length Max. Linear Velocity Void Volume	5% of Total Column Volume 100 Cn 100 Cn/Hr at 45 Pa/g and 4 C 25% Column Volume	Column Equilibration Column Wash Column Regenerate Column Store		Prod. Elution Volume Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	42% Columns Volumes 95% 95% N Y Y Y
T146 Size Excl. Chromatography LPLC	Load Capacity Length Max. Linear Velocity Void Volume	5% of Total Column Volume 100 Cn 100 Cn/Hr at 45 Pa/g and 4 C 25% Column Volume	Column Equilibration Column Wash Column Regenerate Column Store		Prod. Elution Volume Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	42% Columns Volumes 95% 95% N Y Y Y
T147 Dilution	Dilution Factor	3 Liters/Liter	Dilution Time Additional Mix Time		Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	85% 95% Y Y Y Y
T148 Resubstitution	Reagent/Product Ratio Dissolution Time Additional Mix Time	0 L/kg Product 0.50 Hours 0.50 Hours	Reagent 1 Concentration Water Dist.		Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	95% 95% Y Y Y Y

54

Master Process Parameters Table - Biopharmaceutical

Unit Operation Type	Group 1			Group 2			Group 3		
	Parameter	Soln.		Parameter	Soln.		Parameter	Soln.	
749 Enzymatic Modification	Enzyme to Product Ratio Enzyme Concentration Reaction Temp. Reaction Duration	0.084 Liters of Enzyme Stock Per Liter of Start. Proc. Vol. 2 M/Ml 37 Degrees C 30 Minutes 100%		Titration Solution-1 Titration Solution-2 Neutralization	0.007 UL Process 0.02 UL Process 0.57 UL Process		Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	95% 95% Y Y Y	
750 Lyophilization	Product Capacity/Load Product Unit Size	8 Units 100 Grams/Unit		Lyophilization Time Product Weight Reduction	18 Hours 0.85		Step Recovery of Product Step Recovery of T.P. CIP SIP	95% 95% Y Y Y	
751 Heat Exchange	Process Initial Temp. Process Final Temp Utility Initial Temp Utility Final Temp. Process Specific Heat Design Type (P,T,C)	98.6 Degrees C 39.2 Degrees C 34 Degrees C 5 Degrees C 38.6 K BTU/Hr P		Exposure Time	1 Hours		Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	100% 100% Y Y Y	
752 Storage							Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	85% 95% Y Y Y	
753 Fermentation Seed	Scale Up Ratio Fermentor Working Volume Antifoam A Antifoam B Base Acid	10 Fold 50 Liters 1 M/L 1 M/L 5 M/L 5 M/L		Growth Temperature Aeration Sparg Rate Back Pressure Total Duration	37 Hours 1 HP/100L 1.5 YVM 5 PSIG 21 Hrs		Final OD CIP	12 Y	
54 Initial Seeding	Flask Feed Volume Spinner Split Ratio uCarrier Density Number of PBS Washes Number of Media Washes No. of Media/Serum Washes	12 Liters 4 5 Gm/Liter 2 1 2 FBS		Serum Content Feed Rate Days to Confluence	2% FBS 1 Feed per vessel per 2 Days 2 Days		Amplification Factor	1	
55 Culture Vessel Split	Flask Feed Volume Spinner Split Ratio uCarrier Density Number of PBS Washes Number of Media Washes No. of Media/Serum Washes	12 Liters 4 5 Gm/Liter 2 1 2 FBS		Serum Content Feed Rate Days to Confluence	2% FBS 1 Feed per vessel per 2 Days 2 Days		Amplification Factor	1	
56 Culture Flask Split									
57 Stirred Tank Reactor							Step Recovery of Product Step Recovery of T.P. CIP SIP	0.85 95% Y Y Y	
58 Fluidized Bed Reactor	Process Initial Temp. Process Final Temp Utility Initial Temp	37 Degrees C 4 Degrees C 2 Degrees C		Exposure Time	50% Hours		Step Recovery of Product Step Recovery of T.P.	0.85 100%	

A6

Unit Operation Type	Group 1			Group 2			Group 3		
	Parameter	Soln.		Parameter	Soln.		Parameter	Soln.	
	Utility Final Temp. Process Specific Heat Design Type (P,T,C) P		5 Degrees C 12 K BTU/hr P				Temperature Regulation CIP SIP	Y Y Y	
59 Liquid/Liquid Extraction	Liquid/Liquid Ratio Extraction Temperature Addition Duration Additional Mix Duration Mix Energy		1 L Extraction/L Product 4 C 0.5 Hours 4 Hours 0.3 HP/100 L	Phase Separation Time Product Phase (Top/Bottom) Harvest Time		1800% Hours Top 0.5 Hours	Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	Y Y Y Y Y	0.9 50%
60 Solid/Liquid Extraction	Liquid/Liquid Ratio Extraction Temperature Duration Mix Energy		1 L Extraction/L Product 4 C 4 Hours 0.3 HP/100 L	Phase Separation Time Product Phase (Top/Bottom) Harvest Time		1800% Hours Top 0.5 Hours	Step Recovery of Product Step Recovery of T.P. Temperature Regulation CIP SIP	Y Y Y Y Y	0.9 50%